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**'PRODUCT NAME DONOR TREPHINE VACUUM PUNCH**  
**'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

GB - This Instructions for use relates to the following products:  
 DE - Diese Gebrauchsanweisung bezieht sich auf die folgenden Produkte:  
 ES - Estas instrucciones de uso hacen referencia a los siguientes productos:  
 FR - Mode d'emploi des produits suivants :  
 IT - Le seguenti istruzioni per l'uso si riferiscono ai seguenti prodotti:  
 PT - Estas instruções referem-se aos seguintes produtos:  
 RU - Настоящая инструкция по эксплуатации содержит информацию о следующих продуктах:  
 TR - Bu kullanım kılavuzu aşağıda belirtilen ürünler ile ilgilidir:  
 ZH - 本使用说明书涉及以下产品:

GB	Donor Trephine Vacuum Punch
DE	Spenderhornhaut-Trepan-Vakuumstanze
ES	Punzón de succión con trépano para el donante
FR	Emporte-pièce à vide à trépanation donneur
IT	Trapano punch corneale a suzione donatore
PT	Punção a Vácuo para Trépano do Doador
RU	Вакуумный пробойник для донорского трепана
TR	Vakum Punch Donör Trepan
ZH	供体环钻真空冲头

**GB – INSTRUCTIONS FOR USE**

**Description.**

The Donor Trephine Vacuum Punch is a single use only punch and trephine system for the trephination of the donor cornea button. The system consists of a vacuum base, a trephine guide and a standard Coronet trephine (see Fig. 1). It also contains a protector cover as an accessory.

**Intended Use.**

The Donor Trephine Vacuum Punch components allow the surgeon to use a standard trephine to prepare and cut the donor button prior to transplantation. The graduated bowl of the base holds the donor cornea in an anatomically correct position while the trephine guide holds the trephine completely vertical during the procedure so creating a very clean vertical cut through the tissue.

**CAUTION**

- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The punch is supplied sterile and ready for use. Do not use if the pouch has been opened or damaged.
- The Coronet trephine is much sharper than a traditional trephine. Caution should be exercised at all times when handling and using the trephine.
- The Coronet trephine is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface, check the blade integrity under the microscope before use and discard the blade if damaged.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
- This product contains nickel which may cause allergic reaction in patients with sensitivity to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
- For safe disposal, the device should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.
- It is recommended that the trephine be left in the packaging blister until required, to ensure that the blade remains in perfect condition.

**Incident Reporting**

- Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

**STERILISATION**

- The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).
- DO NOT RE-STERILISE. Network Medical Products Ltd. assumes no liability for devices that have been re-sterilised by healthcare facilities.

**Surgical Technique.**

1. Remove the punch base and syringe from the packaging.
2. Remove the cornea from the storage medium by grasping the rim of the surrounding sclera with toothed forceps. Dry the anterior surface with a lint free PVA sponge spear.
3. In preparation for the procedure, it is suggested that an assistant depresses the plunger of the syringe prior to the placement of the donor tissue in the bowl of the vacuum base.
4. Position the corneal tissue in the bowl of the base.
5. Confirm the position of the corneal tissue and if necessary use a micro hook to manipulate and position the tissue optimally.
6. Release the plunger on the syringe to create the vacuum.
7. Once the vacuum has been established, position the trephine guide over the bowl so that it will drop vertically onto the base. Press down gently until the two parts are aligned correctly.
8. Remove the trephine from the packaging. Position the trephine on the guide and gently lower the blade onto the endothelial surface.
9. Press down on the trephine with gentle but firm pressure to push the blade through the corneal tissue until the blade comes to a complete stop. Rotate the trephine carefully through 10 degrees to ensure that the cut is complete.
10. Carefully remove the trephine and the trephine guide together as one piece. The button will remain in the bowl, held gently in place by the vacuum.
11. In the unlikely event that the corneal button remains in the lumen of the trephine, a balanced salt solution or a viscoelastic substance should be gently dripped into the centre of the blade to dislodge the button.
12. Slightly depress the plunger of the syringe or cut the silicone tube to break the vacuum.
13. Cover the corneal tissue with the protector cap.
14. When the recipient cornea has been prepared, the protector cap may be removed and the corneal button removed with a corneal spoon or micro forceps.

**Disposal**

For safe disposal, the trephine should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital protocol.

**HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:**

1. Single use devices have not been validated for re-use. If you re-use a device you may be held **Legally Liable** for the safe performance.
2. Cross-contamination and infection risks to patients. Including transmission of:
  - CJD & Variant CJD.
  - Prion Diseases.
  - Bacterial Endotoxins.
  - Hepatitis B & Hepatitis C.
  - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
  - **Plastics:** Can be weakened, warped or become brittle.
  - **Metals:** Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
  - **Adhesives:** Can degrade if re-sterilised causing parts to dislodge.
  - **Other materials:** May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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**'PRODUCT NAME DONOR TREPHINE VACUUM PUNCH  
'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

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