LONG HANDLED TREPINE

PRODUCT INFORMATION AND INSTRUCTIONS

 **NETWORK MEDICAL PRODUCTS LTD**

Coronet House, Kearsley Road,
Ripon, North Yorkshire, HG4 2SG, UK

Tel: +44 (0)1765 609555 Fax: +44 (0)1765 608476

info@networkmedical.co.uk | www.networkmedical.co.uk

 **1639**

GB - This Instructions for use relates to the following products:

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GB - Instructions for Use

Description

The Coronet trephine is a single use only device with an ultra-sharp cutting edge, supplied with a knurled finger grip for secure manual trephination.

Intended Use

The CORONET range of corneal trephines are devices used in Corneal Graft Surgery and can be used for full penetrating keratoplasty of the recipient cornea or for the preparation of the donor button during corneal transplant procedures. The ultra-sharp, thin profile cutting edge minimises damage to the endothelial cells of the cornea and reduces undercut and tissue distortion.

CAUTION

- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- This device is supplied sterile and ready for use. Do not use if the pouch has been opened or damaged.
- The Coronet trephine is much sharper than a traditional trephine. Caution should be exercised at all times when handling and using the trephine.
- The Coronet trephine is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface, check the blade integrity under the microscope before use and discard the blade if damaged.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
- This product contains nickel which may cause allergic reaction in patients with sensitivity to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
- For safe disposal, the device should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.
- To avoid damaging the blade edge, take care when removing the device from the rigid plastic container. Using finger and thumb grasp the handle between the two holding mounts and lift the trephine vertically.

Incident Reporting

- Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

STERILISATION

- The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).
- **DO NOT RE-STERILISE.** Network Medical Products Ltd. assumes no liability for devices that have been re-sterilised by healthcare facilities.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

1. Single use devices have not been validated for re-use. If you re-use a device you may be held **Legally Liable** for the safe performance.
2. Cross-contamination and infection risks to patients. Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
 - **Metals:** Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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Advena Limited