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



GB - This Instructions for use relates to the following products:

PRODUCT NAME ENDOGLIDE ULTRATHIN
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

DE - Diese Gebrauchsanweisung bezieht sich auf die folgenden Produkte:
 ES - Estas instrucciones de uso hacen referencia a los siguientes productos:
 FR - Mode d'emploi des produits suivants :
 IT - Le seguenti istruzioni per l'uso si riferiscono ai seguenti prodotti:
 NL - Deze gebruiksaanwijzing heeft betrekking op de volgende producten:
 PT - Estas instruções referem-se aos seguintes produtos:
 TR - Bu kullanım kılavuzu aşağıda belirtilen ürünler ile ilgilidir:

GB, DE, ES, FR, IT, NL, PT, TR	Endoglide Ultrathin
--------------------------------	---------------------

GB – INSTRUCTIONS FOR USE

Description

The EndoGlide Ultrathin (UT) is a device designed particularly to meet the needs of the corneal graft surgeon performing DSEK (Descemet's Stripping Endothelial Keratoplasty) or DSAEK (Descemet's Stripping Automated Endothelial Keratoplasty) procedures using ultrathin donor tissue. The device is a sterile four-part single use only delivery system for the previously prepared donor tissue, which allows the surgeon to hold and deliver the tissue in a manner which minimises any endothelial cell damage and consequent loss of function.

Intended Use

The device is intended solely for the delivery and insertion of previously prepared donor cornea tissue for transplantation during a DSEK or DSAEK procedure.

CAUTIONS

- This device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the blister pack has been opened or damaged.
- Place the preparation base on a flat, sterile surface before use.
- Carefully place the donor cornea, comprising the posterior lenticule only, endothelial surface up within the well of the Preparation Base before proceeding further.
- Before removing the cartridge from the preparation base ensure that the locking mechanism has engaged correctly and that the two parts are securely locked together.
- CAUTION should be exercised when inserting the cartridge if the patient has an intraocular lens implant.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.

Possible Adverse Effects of the Surgical Procedure

- Complications which may occur include but are not limited to: Graft detachment/rejection/ failure, Glaucoma, Retinal detachment, macular hole, epiretinal membrane, raised intraocular pressure, corneal synechia and pupillary block. Donor-host interface vascularization may occur.
- In some phakic eyes, mild anterior crystalline lens opacities have also been reported. Subepithelial haze and interface pigment deposits have also been reported
- Small areas of interface fluid at the edge of the graft have been known to occur with DSEK/DSAEK procedures.
- With very thin tissue, care should be taken as folds may occur. Partial dehiscence of the posterior lamellar graft may occur.

Incident Reporting

- Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Protocol for surgical technique

Recipient cornea preparation

1 Wound construction

A temporal scleral or clear corneal incision may be used to insert the EndoGlide UT containing the donor tissue. For scleral incisions, a 4.5 mm blade should be used. For clear corneal incisions, a 4.9mm blade should be used. In both cases, the wound depth should be at least 1.5mm deep to create a self-sealing wound.

2 AC maintainer

Use an AC maintainer to ensure formation of the AC during insertion. Placement of the AC maintainer could be adjacent to either the superior or inferior aspect of the temporal incision so as not to obstruct the donor incision. BSS flow through the maintainer should be away from the incision and low to moderate during the entire insertion procedure.

3 Nasal paracentesis

Perform a 1mm paracentesis at the nasal limbus or clear cornea opposite to the temporal incision through which to introduce the Tan EndoGlide Placement Forceps (TEP Forceps) to pull the donor tissue into the AC.

Donor Cornea Preparation and Sizing

- 1 The EndoGlide UT is specifically designed for loading and insertion of an ultrathin posterior lenticule prepared by a microkeratome. (Note : the EndoGlide UT can be used with lenticules of thickness 70 to 250µm).
- 2 Donor trephination with a suitable Coronet Donor Punch: the Coronet EndoGlide UT is designed to enable insertion of a corneal donor up to 9.5mm in diameter.

Cartridge preparation

The procedure should be performed under the operating microscope with the preparation base positioned on a flat, stable, sterile surface. Prepare the recipient cornea as above, so that the donor tissue insertion may be performed immediately after the donor is coiled in the cartridge. Note that the cartridge is supplied pre-loaded onto the preparation base.

1 Positioning of EndoGlide Cartridge

The cartridge and saddle are supplied pre-loaded on the preparation base.

2 Donor tissue placement

Maximal forward placement of the donor to abut the cartridge opening makes it easier to grasp the donor for the coiling process. This brings the donor edge closer to the loading forceps, and positions the donor at the ideal height for easy grasping with the forceps for the pull-through. The forceps guide bridge should be used to aid the correct placement of the forceps.

There are two methods for tissue placement – a. for use with a single lenticule, or b. for use when supplied with both the anterior and posterior lenticules.

a. Single Lenticular Method

- Place a small amount of BSS in the well of the preparation base for lubrication. Inject BSS into the inner chamber of the cartridge, which can be completely filled. Using a Paton spatula, carefully place the donor cornea endothelial surface up in the well of the preparation base. Apply the BSS cannula to the edge of the lenticule, and gently push the lenticule close to the entrance of the saddle.
- After correct positioning of the donor, place a very thin stream of dispersive viscoelastic onto the donor surface, along the long axis of the cartridge, prior to donor coiling. This provides for additional endothelial protection during the pull-through process.

b. Double (anterior and posterior) Lenticular Method

- Place a small amount of BSS in the well of the preparation base and carefully inject BSS into the inner chamber of the cartridge, which should be completely filled. Using a Paton spatula, carefully place the donor cornea, comprising both the anterior and posterior lenticules, endothelial surface up. Physical separation may be gently performed via gentle hydro-dissection of the two lenticules. Apply the BSS cannula to the edge of the lenticule, and gently push the lenticule close to the entrance of the saddle.
- After correct positioning of the donor, place a very thin stream of dispersive viscoelastic onto the donor surface, along the long axis of the cartridge, prior to donor coiling. This provides for additional endothelial protection during the pull-through process.

3 Curling procedure

Introduce the Tan EndoGlide Loading Forceps (TEL forceps) through the anterior opening of the cartridge into the chamber up to the level of saddle, where the leading edge of the donor will be positioned. The forceps guide bridge should be used to aid the correct placement of the forceps. Using the TEL forceps grasp the edge of the donor, and gradually pull it into the cartridge chamber. As the donor tissue enters the cartridge, the sides of the donor tissue will naturally curl upwards, without wrinkling, along the upward curving walls of the saddle tunnel, and the tissue will form a double coil configuration in the cartridge. It is usually unnecessary to use a cannula to stroke the edges of the donor upwards along the side walls of the chamber. Pull the donor tissue all the way in, until the leading anterior edge just reaches the anterior opening of the cartridge. Gently release the donor tissue from the forceps and remove the forceps out of the cartridge. If the tissue remains attached to the teeth of the forceps, slight rotational movements of the forceps along its long axis will facilitate tissue release. The donor should now be fully formed in a double coil configuration within the cartridge. **NOTE:** This is inferior at this stage as the cartridge is positioned upside down on the preparation base.

The EndoGlide Introducer

1 Removal of Donor Anterior Lenticule (Double Lenticular Method only)

Remove the donor anterior lenticule from the well of the Preparation Base and discard.

2 Removal of the Saddle

Holding the 'saddle' between finger and thumb, release the component one side at a time and then lift off vertically to remove the saddle before positioning the introducer.

3 Attachment of the EndoGlide Introducer

PRODUCT NAME ENDOGLIDE ULTRATHIN
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

The introducer slides into position behind the cartridge using the guiding grooves. Advance the introducer, with word "loading" facing uppermost, into the posterior end of the cartridge until it fully engages with an audible click.

4 Removal of the EndoGlide Cartridge, Introducer and donor complex

The loaded cartridge, with the introducer attached, may now be removed in its' entirety from the preparation base by sliding the assembly backwards and out of the grooves of the base. The cartridge is now ready for insertion into the recipient eye.

EndoGlide Cartridge Insertion

NOTE: After removing the introducer and cartridge from the base, turn the assembly right-side up (the word "insertion" printed on the introducer should now be uppermost) to ensure correct orientation for insertion into the recipient eye.

- 1 Grasp the introducer right way up with the thumb and forefinger for easy insertion into the wound (thumb uppermost) (similar to insertion of a thumbdrive into the USB port of a computer).
- 2 With the AC maintainer on low flow, gently slide the anterior tongue (glide) of the cartridge through the temporal wound, to advance the cartridge into the AC. The flat surface of the tongue ensures moderate closure of the wound to reduce egress of AC fluid and prevents iris prolapse through the wound.
- 3 Grip the anterior lip of the wound (to prevent lip inversion), and introduce the sloping anterior opening of the cartridge fully into the wound until the uppermost edge of the opening is completely engaged through the wound. Advancement should be fairly rapid to reduce chamber collapse. Continue to advance the cartridge further until the anterior opening is seen through the cornea to be fully within the AC. At this point the forward edge of the cartridge tongue (glide) is positioned across the pupil and over the nasal iris, and the EndoGlide should snugly seal the wound.

Donor Pull-Through and Removal of the EndoGlide

With one hand still holding the introducer and cartridge in position (the EndoGlide can be used to steady the eye at this point), insert the TEP Forceps into the AC through the nasal paracentesis with the other hand. Advance the forceps into the AC, over the cartridge glide, into the anterior opening of the cartridge to reach the leading stromal edge of the donor tissue. Grasp the donor and gently pull the donor tissue out of the cartridge chamber into the AC. Keep the cartridge and introducer in position during this stage. Whilst still holding onto the donor tissue with the TEP Forceps, retract the cartridge and introducer out of the eye. The wound will close, and the AC will remain deep due to the AC maintainer. Once the donor tissue is fully placed in the AC, it will automatically start to uncurl in the anatomically correct endothelial down position. Gentle sideways "to-and-fro" shaking of the donor with the forceps will easily assist in uncurling if one wing remains slightly curled. During this process, the AC should remain deep and fully formed with the AC maintainer still on. If the edges of the donor are caught on the iris, and do not fully uncoil, deepening the chamber by increasing AC maintainer flow may help achieve full uncurling. Gently tapping the cornea from the corneal surface may also assist in full donor uncurling. **Note: do not release the donor with the forceps at this stage.**

Donor centration and air injection

- 1 Ensure that the donor is centrally placed over the recipient cornea with the forceps. Inject a small amount of air (2mm size) with an air cannula beneath the donor – the air bubble ensures that the donor continues to float up against the recipient stromal surface. Do not attempt to insert a large air bubble at this stage. If the donor is not completely central, a large bubble will quickly press the donor up against the recipient stromal surface, initiating sticking and make subsequent donor positioning difficult.
- 2 Finally, release the donor from the TEP forceps and retract the forceps out of the eye. If there is concern that the wound is not fully self-sealing, and there is high vitreous pressure (in which case the AC may shallow and the donor may shoot out through the wound if it is released) then before releasing the donor, the AC maintainer may be shut off or reduced, and the lip of the wound may be held down with forceps before releasing the donor fully. Full wound closure, air tamponade and completion of the DSAEK procedure should then be performed in the usual manner.

Sterilisation

- The EndoGlide Ultrathin is supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO)

Hazards associated with RE-USE of this Single Use Only device:

1. This single use device has not been validated for re-use. If you re-use a device you may be held **Legally Liable** for the safe performance.
2. Cross-contamination and infection risks to patients, including the transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS

Note : Reprocessing may not remove all viable micro-organisms. Abnormal proteins associated with prion diseases eg. Creutzfeldt-Jakob disease (CJD) and variant Creutzfeldt-Jakob disease, have been identified in corneal tissue and are very resistant to all conventional methods of decontamination.
3. Device failure through material fatigue or degradation caused by initial use and design. Plastics can be weakened, warped or become brittle. The accuracy, function and performance of the device **will be seriously affected if the device is re-sterilised.**
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.
5. **Separating components** - Forcibly separating the introducer from the cartridge after primary use will lead to leakage of aqueous from the back of the cartridge and loss of intraocular pressure in any re-use.

Also Recommended for use with the Endoglide Ultrathin

53-951	Tan EndoGlide Placement Forceps
53-952	Tan EndoGlide Loading Forceps

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