**CRAWFORD
INTUBATION SET &
CANALICULUS
INTUBATION SET**

**PRODUCT
INFORMATION AND
INSTRUCTIONS**

NETWORK MEDICAL PRODUCTS LTD

Coronet House, Kearsley Road,
Ripon, North Yorkshire, HG4 2SG, UK

Tel: +44 (0)1765 609555 Fax: +44 (0)1765 608476

info@networkmedical.co.uk | www.networkmedical.co.uk

GB - This Instructions for use relates to the following products:

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ES - This instructions for use relates to the following products:

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
GB - Instructions for Use

Description

The **NETWORK MEDICAL Crawford Intubation Set and Canaliculus Intubation Set** consist of two malleable stainless steel probes attached by a length of silicone tubing. The probes of the Crawford Intubation Set have olive shaped tips which have been designed to prevent unintentional perforation of the soft tissue and can be located by the hook of the Crawford Retrieval device and used to draw the probe down through the nose efficiently and safely. Both sets are supplied sterile for **SINGLE USE ONLY**.

Intended Use

The **NETWORK MEDICAL Crawford Intubation and Canaliculus Intubation Sets** are intended for silicone intubation procedures. The device provides a means of restoring the patency of the lacrimal system.

CAUTIONS

- This device is supplied **STERILE** and ready to use.
- This device is for **SINGLE USE ONLY**. Do **NOT** re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is only to be used by a suitably trained or qualified person under the facility approved and recommended conditions.
- Ease of passage of the probes through the Lachrymal System will vary widely from patient to patient with the placement of the tubes in some patients being quite difficult due to narrow openings in the lachrymal system, especially at the lower end of the bony canal. Separation of the silicone tubing from the probe may occur when intubating patients with severe blockage or a pronounced bony prominence.
- Silicone tubing could be difficult to detect or to remove. Examination under anaesthesia, including nasal endoscopy may be warranted when accidental tube extrusion cannot be reliably confirmed.
- For safe disposal, the needle should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the silicone tubing after use to prevent contamination.
- Do not use for any purposes other than lacrimal surgery.

MRI Compatibility

All tubing is made from Silicone, which is MRI compatible.

Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Instructions and Surgical Procedure

- Each punctum should be dilated.
- The lachrymal system is probed to open up any blockages using a standard lachrymal probe such as a Bowman 00.
- The Crawford probe is passed through the Upper Punctum and across the Upper

Canaliculus and then orientated down through the lachrymal system into the nose approximately 35mm – 40mm.

- If any difficulties are experienced in locating the probe a Bowman probe may be used to establish metal to metal contact and define the position of the probe. The surgeon must keep lateral tension on the lid to avoid an accordion type effect of the canaliculus, which increases the risk of perforating the epithelium and creating a false passage.
- The probe can be located in the nose, lateral to the Inferior Turbinate in the Inferior Meatus of the nose using the Crawford Retrieval Device. The hook of the Retrieval Device is located over the olive tip of the probe and the probe pulled from the nose. Often a push / pull technique will be required to complete the intubation.
- If using the Canaliculus Intubation set, the probes are located and pulled from the nose using forceps. Often a push / pull technique will be required to complete the intubation.
- The second probe is passed down the Inferior Canaliculus and out of the nose in a similar fashion.
- With the silicone tubing in place the probes are removed and tubing tied and the knotted tubing secured by a fine suture of the surgeon's choice. The knotted tubes are then tucked up into the nose.
- The device is made from implantable grade silicone and the tubing is typically removed from 2 weeks to 15 months as deemed appropriate by the surgeon.
- Removal of the tubes is accomplished by locating the tube in the eye between the puncta and pulling the tube upwards. The tubing can then be cut and removed.
- The tubes must be removed through the nose. In older patients, it may be possible to expel the knot by having the patient sneeze, but in small children it may require a second general anaesthetic to locate the knot and remove the tubes.

Possible Adverse Effects

- The device is made from implantable grade silicone however the surgeon should regularly monitor the patient for suitability of the implant. The implant should be removed or replaced if damaged when deemed appropriate by the surgeon.
- Complications which may occur include, but are not limited to: granuloma formation, cellulitis, infection, epistaxis, punctal or canalicular erosion, corneal erosion or abrasion, asymptomatic cheesewiring or marsupialization, nasal discharge, periorbital edema and periorbital ecchymosis.
- It is important to note that Canaliculi or puncta may become lacerated or split if the tubes are too tight. Care must be taken in children to allow for patient growth. The surgeon should regularly monitor the patient for any signs of damage around the implant for punctal erosion, elongation or any sign of tube migration through the puncta.
- Early post-surgical displacement of the implant may occur particularly in children or when forcibly sneezing, potentially leading to lid notching or development of an ectropion due to poor healing.
- Patients requiring Continuous Positive Airway Pressure (CPAP) therapy may experience irritation following the procedure due to air regurgitation.
- Particular attention should be taken during and prior to knotting/suturing to avoid over extension of tubing.
- There is a risk of synechia formation between the lateral nasal wall and the middle turbinate mucosa.
- Tubing is not X-Ray or CT Scan detectable.
- The probe contains nickel which may cause allergic reaction in patients with sensitivity

to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.

- This product contains silicone which may cause allergic reaction in patients with sensitivity to silicone. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to silicone.
- There is a potential risk of synechia formation between the lateral nasal wall and the middle turbinate mucosa
- There is a potential risk of cross infection if devices are not disposed of after use following hospital approved procedures.
- There is a risk of periorbital edema and/or periorbital ecchymosis.

Sterilisation

The NETWORK MEDICAL Crawford Intubation Set, Retrieval Device and Canaliculus Intubation sets are supplied sterile by Gamma Irradiation and are designated for SINGLE USE ONLY.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES

1. Single use devices have not been validated for re-use.
 - If you re-use a device you may be held Legally Liable for the safe performance.
2. Cross-contamination and infection risks to patients . Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
 - Plastics: Can be weakened, warped or become brittle.
 - Metals: Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
 - Adhesives: Can degrade if re-sterilised causing parts to dislodge.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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Advena Limited