
'PRODUCT NAME DALK CLEAR-VUE™ RECIPIENT VACUUM TREPINE
'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

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'PRODUCT NAME DALK CLEAR-VUE™ RECIPIENT VACUUM TREPHINE
'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

GB - This Instructions for use relates to the following products:
 DE - Diese Gebrauchsanweisung bezieht sich auf die folgenden Produkte:
 ES - Estas instrucciones de uso hacen referencia a los siguientes productos:
 FR - Mode d'emploi des produits suivants :
 IT - Le seguenti istruzioni per l'uso si riferiscono ai seguenti prodotti:
 PT - Estas instruções referem-se aos seguintes produtos:
 RU - Настоящая инструкция по эксплуатации содержит информацию о следующих продуктах:
 TR – Bu kullanım kılavuzu aşağıda belirtilen ürünler ile ilgilidir:
 ZH - 本使用说明书涉及以下产品:

GB	DALK Clear-Vue™ Recipient Vacuum Trephine
DE	DALK Clear-Vue™ Empfängerhornhaut-Vakuum-Trepan
ES	DALK Clear-Vue™ Recipient Vacuum Trephine
FR	Trépan à vide receveur Clear-Vue™ DALK
IT	Trapano corneale a suzione ricevente DALK Clear-Vue™
PT	Trépano de Vácuo para Recetor DALK Clear-Vue™
RU	Вакуумный трепан для реципиента DALK Clear-Vue™
TR	DALK Clear-Vue™ Vakumlu Alici Trefini
ZH	DALK Clear-Vue™ 受体真空环钻

GB – INSTRUCTIONS FOR USE

Description

The CORONET® DALK Clear-Vue™ Recipient Vacuum Trephine has an ultra-sharp circular blade which is placed inside the DALK Body. The centering tube is supplied already housed within the DALK body and attached to a spring-loaded syringe. Evacuating the syringe when the trephine is placed in contact with the cornea creates a vacuum which holds the trephine firmly in place on the cornea. Crosshairs aid centration over the optical axis.

Note

The centering tube is supplied already housed within the DALK body. Parallax is minimised as the x-hairs are situated at the very distal end of the centering tube

Intended Use

The CORONET range of corneal trephines are devices used in Corneal Graft Surgery and can be used for full penetrating keratoplasty of the recipient cornea or for the preparation of the donor button during corneal transplant procedures. The ultra-sharp, thin profile cutting edge minimises damage to the endothelial cells of the cornea and reduces undercut and tissue distortion. The DALK Clear-Vue Recipient Vacuum Trephine is a single use trephine designed for DALK (Deep anterior lamellar keratoplasty).

CAUTION

- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- This device is supplied sterile and ready for use. Do not use if the pouch has been opened or damaged.
- The Coronet trephine is much sharper than a traditional trephine. Caution should be exercised at all times when handling and using the trephine.
- The Coronet trephine is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface, check the blade integrity under the microscope before use and discard the blade if damaged.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician
- This product contains nickel which may cause allergic reaction in patients with sensitivity to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
- For safe disposal, the device should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.
- It is recommended that the trephine be left in the packaging blister until required, to ensure that the blade remains in perfect condition.

Incident Reporting

- Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Surgical Technique

- Evacuate the syringe and hold the syringe in situ.
- Place the DALK body into position and place in the desired location through using the centering tube. The body is pressed into position to assist the formation of the vacuum.
- Release the syringe
- Remove the Centering tube and discard
- Carefully insert the DALK guarded trephine into the DALK body and gently advance onto the surface of the cornea
- Using gentle pressure, advance the guarded trephine until resistance is felt.
- Cut the silicone tube to break the vacuum and remove the DALK body assembly from the recipient cornea.
- Continue the procedure as intended.

STERILISATION

- The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).
- DO NOT RE-STERILISE.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

- Single use devices have not been validated for re-use. If you re-use a device you may be held **Legally Liable** for the safe performance.
- Cross-contamination and infection risks to patients. Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
- Device failure through material fatigue or degradation caused by initial use and design:
 - Plastics: Can be weakened, warped or become brittle.
 - Metals: Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
 - Adhesives: Can degrade if re-sterilised causing parts to dislodge.
 - Other materials: May degrade, becoming unacceptable when compared to original manufacturing criteria.
- Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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