**PRODUCT NAME DMEK GUARDED TREPHINE PUNCH**  
**PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

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

**PRODUCT NAME DMEK GUARDED TREPINE PUNCH**  
**PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

GB - This Instructions for use relates to the following products:  
 DE - Diese Gebrauchsanweisung bezieht sich auf die folgenden Produkte:  
 ES - Estas instrucciones de uso hacen referencia a los siguientes productos:  
 FR - Mode d'emploi des produits suivants :  
 IT - Le seguenti istruzioni per l'uso si riferiscono ai seguenti prodotti:  
 PT - Estas instruções referem-se aos seguintes produtos:  
 RU - Настоящая инструкция по эксплуатации содержит информацию о следующих продуктах:  
 TR - Bu kullanım kılavuzu aşağıda belirtilen ürünler ile ilgilidir:  
 ZH - 本使用说明书涉及以下产品:

GB	DMEK Guarded Trephine Punch
DE	Geschlossene DMEK-Trepanstanz
ES	Punzón trepanador protegido para DMEK
FR	Punch Trépan protégé DMEK
IT	Trapano punch corneale con protezione DMEK
PT	Punção para Trépano DMEK com Capa
RU	Трепан-пробойник с защитой для DMEK
TR	DMEK Korumalı Trepan Punch
ZH	DMEK 防护式环钻冲头

**GB – INSTRUCTIONS FOR USE**

**Description**

The DMEK Guarded Trephine Punch is a single use only punch and trephine system for the trephination of the donor cornea for DMEK. The system consists of a vacuum base, a trephine guide and a guarded CORONET™ trephine (see Fig. 1).

**Intended Use**

The CORONET range of corneal trephines are devices used in Corneal Graft Surgery and can be used for full penetrating keratoplasty of the recipient cornea or for the preparation of the donor button during corneal transplant procedures. The ultra-sharp, thin profile cutting edge minimises damage to the endothelial cells of the cornea and reduces undercut and tissue distortion. The DMEK Guarded Trephine Punch is designed to allow the donor corneal tissue to be held in the most relaxed state in the bowl, while a vacuum is drawn. The vacuum holds the tissue gently in the optimal position, in preparation for the cutting of the donor cornea. The guarded trephine, supplied as part of the system, is ultra-sharp and makes a clean vertical cut into the corneal tissue.

**CAUTION**

- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- This device is supplied sterile and ready for use. Do not use if the pouch has been opened or damaged.
- The Coronet trephine is much sharper than a traditional trephine. Caution should be exercised at all times when handling and using the trephine.
- The Coronet trephine is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface, check the blade integrity under the microscope before use and discard the blade if damaged.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
- This product contains nickel which may cause allergic reaction in patients with sensitivity to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
- For safe disposal, the device should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.
- It is recommended that the trephine be left in the packaging blister until required, to ensure that the blade remains in perfect condition.
- Partial penetration of the donor button requires minimal pressure on the trephine.

**Incident Reporting**

- Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

**Surgical Technique**

1. Remove the punch base and syringe from the packaging.
2. Remove the cornea from the storage medium by grasping the rim of the surrounding sclera with toothed forceps. Dry the anterior surface with a lint free PVA sponge spear.
3. In preparation for the procedure, it is suggested that an assistant depresses the plunger of the syringe prior to the placement of the donor tissue in the bowl of the vacuum base.
4. Position the corneal tissue in the bowl of the base.
5. Confirm the position of the corneal tissue and if necessary use a micro hook to manipulate and position the tissue optimally.
6. Release the plunger on the syringe to create the vacuum.
7. Once the vacuum has been established, position the trephine guide over the bowl so that it will drop vertically onto the base. Press down gently until the two parts are aligned correctly.
8. Remove the guarded trephine from the packaging. Position the trephine in the guide and gently lower the blade onto the surface of the endothelium.
9. Gently tap on all four quadrants of the trephine.

**Or**

- 9B. Place finger or thumb over the top of the trephine to ensure equal pressure. Press firmly until resistance is felt. These steps are required to achieve a 360° cut.
10. Carefully withdraw the trephine and remove from the guide.
11. Carefully detach the Trephine Guide from the base and remove.
12. The Donor Cornea is now ready for final preparation of the DMEK lenticule, in the desired manner.
13. Finally, when required, slightly depress the syringe or cut the silicone tube to break the vacuum.

**STERILISATION**

- The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).
- DO NOT RE-STERILISE. Network Medical Products Ltd. assumes no liability for devices that have been re-sterilised by healthcare facilities.

**HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:**

1. Single use devices have not been validated for re-use. If you re-use a device you may be held Legally Liable for the safe performance.
2. Cross-contamination and infection risks to patients. Including transmission of:
  - CJD & Variant CJD.
  - Prion Diseases.
  - Bacterial Endotoxins.
  - Hepatitis B & Hepatitis C.
  - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
  - Plastics: Can be weakened, warped or become brittle.
  - Metals: Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
  - Adhesives: Can degrade if re-sterilised causing parts to dislodge.
  - Other materials: May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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