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

**PRODUCT NAME DMEK ENDOGLIDE**  
**PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

|    |           |                                |              |                       |                          |
|----|-----------|--------------------------------|--------------|-----------------------|--------------------------|
| FR | Cartouche | Pont de guidage des forceps    | Introduteur  | Base de préparation   | Puits de préparation     |
| IT | Cartuccia | Ponte guida pinza              | Introduttore | Base di preparazione  | Pozzetto di preparazione |
| PT | Cartucho  | Conector Guia do Fórceps       | Introdutor   | Base de Preparação    | Poço para Preparação     |
| RU | Картридж  | Направляющая перемычка пинцета | Интродьюсер  | Основа для подготовки | Подготовительная лунка   |
| TR | Kartuş    | Forseps kılavuz köprüsü        | Introdüser   | Preparat Kaidesi      | Preparat Haznesi         |

GB - This Instructions for use relates to the following products:  
 DE - Diese Gebrauchsanweisung bezieht sich auf die folgenden Produkte:  
 ES - Estas instrucciones de uso hacen referencia a los siguientes productos:  
 FR - Mode d'emploi des produits suivants :  
 IT - Le seguenti istruzioni per l'uso si riferiscono ai seguenti prodotti:  
 PT - Estas instruções referem-se aos seguintes produtos:  
 RU - Настоящая инструкция по эксплуатации содержит информацию о следующих продуктах:  
 TR - Bu kullanım kılavuz aşağıda belirtilen ürünler ile ilgilidir:

|    |                |
|----|----------------|
| GB | DMEK EndoGlide |
| DE | DMEK EndoGlide |
| ES | DMEK EndoGlide |
| FR | EndoGlide DMEK |
| IT | DMEK EndoGlide |
| PT | EndoGlide DMEK |
| RU | DMEK EndoGlide |
| TR | DMEK EndoGlide |

**GB – INSTRUCTIONS FOR USE**

**Description**

The DMEK EndoGlide is a device designed particularly to meet the needs of the corneal graft surgeon performing DMEK (Descemet Membrane Endothelial Keratoplasty) procedures. The device is a sterile three-part single use only delivery system for the previously prepared donor tissue, which gives the surgeon enhanced donor tissue control and stabilisation of the anterior chamber during donor graft insertion and minimises any endothelial cell damage and consequent loss of function.

**Intended Use**

The device is intended solely for the delivery and insertion of previously prepared donor cornea tissue for transplantation during a DMEK procedure..

**CAUTIONS**

- This device is supplied STERILE and ready to use.
- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device
- After use, the device should be disposed of in accordance with hospital protocol.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
- Place the preparation base on a flat, sterile surface before use.
- CAUTION should be exercised when inserting the cartridge if the patient has an intraocular lens implant.
- Place the corneal button with the tri-folded DMEK graft onto the Preparation well and orientate the trifold graft edge to align with the cartridge opening.
- Before removing the cartridge from the preparation base, ensure that the introducer has engaged correctly and that the two parts are securely locked together.
- The Use of Viscoelastic is not recommended.

**Possible Adverse Effects of the Surgical Procedure**

- Complications which may occur include but are not limited to: Graft detachment/rejection/ failure, Glaucoma, Retinal detachment, macular hole, epiretinal membrane, raised intraocular pressure, corneal synechia and pupillary block. Donor-host interface vascularization may occur.
- In some phakic eyes, mild anterior crystalline lens opacities have also been reported. Subepithelial haze and interface pigment deposits have also been reported
- Small areas of interface fluid at the edge of the graft have been known to occur with DSEK/DSAEK procedures.
- With very thin tissue, care should be taken as folds may occur. Partial dehiscence of the posterior lamellar graft may occur.
- There is a risk of a mild hyperopic shift with DMEK.

**Incident Reporting**

- Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

**Protocol for surgical technique**

**Recipient cornea preparation**

**1. Wound construction**

A 2.65mm blade is used to create a clear corneal incision through which the EndoGlide cartridge may be introduced. The wound depth should be at least 1.5mm deep to create a self-sealing wound.

**2. Nasal paracentesis**

Perform a 1mm paracentesis at the nasal limbus or clear cornea opposite to the temporal incision through which to introduce the Tan EndoGlide Placement Forceps (TEP forceps) to pull the donor tissue into the AC.

**3. AC maintainer**

Use an AC maintainer to ensure formation of the AC during insertion. Placement of the AC maintainer could be adjacent to either the superior or inferior aspect of the nasal paracentesis incision so as not to obstruct the donor incision. BSS flow through the maintainer should be away from the incision and low to moderate during the entire insertion procedure.

**Donor Cornea Preparation**

- Fully stripped DMEK tissue should be laid flat endothelium up on corneal button
- Perform tri-fold folding with a small amount of BSS fluid present
- Use Kelman forceps to slide the tri-folded donor tissue out towards the scleral edge.

**Cartridge preparation and Sizing**

The procedure should be performed under the operating microscope with the preparation base positioned on a flat, stable, sterile surface. Prepare the recipient cornea as above, so that the donor tissue insertion may be performed immediately after the cartridge is ready. Note that the cartridge is supplied assembled on the preparation base.

**1. Positioning of EndoGlide Cartridge**

The cartridge is supplied assembled on the preparation base.

**2. Donor tissue placement: Tri-fold loading procedure**

- Place the corneal button with tri-folded DMEK graft onto the preparation well. Orientate the trifold graft edge to align with the cartridge opening.
- Ensure that there is BSS between the cartridge opening and the DMEK donor edge
- Use the forceps guide bridge, pass the loading forceps through the cartridge and depress the scleral edge of the donor button to bring the DMEK graft directly into the forceps tips.
- Grasp the DMEK graft with the loading forceps and gently advance the graft into the cartridge (keep sclera depressed until graft is inside the cartridge). The anterior edge of the graft should reach the anterior opening of the cartridge.

**The EndoGlide Introducer**

**1. Removal of donor stromal button**

Remove the donor corneal/stromal button from the well of the preparation base and discard.

**PRODUCT NAME DMEK ENDOGLIDE**  
**PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

- Attachment of the EndoGlide Introducer**  
The introducer slides into position behind the cartridge using the guiding grooves. Advance the introducer, with word "loading" facing uppermost, into the posterior end of the cartridge until it fully engages with an audible click.
- Removal of the EndoGlide Cartridge, Introducer and donor DMEK graft**  
The loaded cartridge, with the introducer attached, may now be removed in its' entirety from the preparation base by sliding the assembly backwards and out of the grooves of the base. The cartridge is now ready for insertion into the recipient eye.

**EndoGlide Cartridge Insertion**

**NOTE:** After removing the introducer and cartridge from the base, turn the assembly right way up (the word "insertion" printed on the Introducer should now be uppermost) to ensure correct orientation for insertion into the recipient eye.

- Grasp the introducer right side up with the thumb and forefinger for easy insertion into the wound (thumb uppermost) – Note this is similar to insertion of a thumbdrive into the USB port of a laptop computer.
- Gently insert the anterior tongue (glide) of the cartridge into the temporal wound, and advance into the anterior chamber (AC). The flat surface of the tongue ensures moderate closure of the wound to reduce egress of AC fluid and prevents iris prolapse through the wound. Moderate flow of BSS from the AC maintainer ensures maintenance of the AC.
- Introduce the sloping anterior opening of the cartridge through the wound and advance it fully until the uppermost edge of the opening is completely engaged through the wound. More rapid advancement at this stage ensures minimal egress of BSS and maintenance of the AC. Advance the cartridge further until the anterior opening is seen through the cornea to be fully within the AC. At this point the forward edge of the cartridge tongue (glide) will be over the nasal iris.

**Donor Pull-Through**

- With one hand still holding the introducer and cartridge in position (the EndoGlide can be used to steady the eye at this point), insert the Tan EndoGlide Placement forceps (TEP Forceps) into the AC through the nasal paracentesis with the other hand. Advance the forceps into the AC, over the cartridge tongue (glide), into the anterior opening of the cartridge to reach the leading stromal edge of the donor tissue.
- Use the forceps to grasp the stromal edge of the donor tissue and gently pull the donor tissue out of the cartridge chamber into the AC. Keep the cartridge and introducer in position during this stage. Once the donor tissue is fully placed in the AC, it will automatically start to uncurl into position, endothelial side down, with both wings of the donor uncurling.
- Gentle sideways "to-and-fro" shaking of the donor with the forceps will easily assist in uncurling if one wing remains slightly curled. During this process, **the ACM flow must be reduced to very low** (or turned off) to ensure a shallow, non-turbulent AC during the donor pull-through.
- Donor pull through should be very slow. Ensure a good grasp of the leading edge of the graft and gradually advance it into the AC. Ensure that the AC remains shallow but stable throughout. Once the donor is ¾ out of the cartridge, and has completely moved beyond the corneal wound, the cartridge can be slowly eased out of the wound (i.e. avoid leaving the back edge of the graft in the corneal wound.)

Note: avoid a very tight wound as this will make it difficult to remove the cartridge without pulling the eye, which may then detach the ACM. If this occurs the ACM must be reattached.

**Removal of the EndoGlide, donor centration and air injection**

- Whilst still holding onto the donor graft tissue with the TEP Forceps, retract the cartridge and introducer assembly out of the eye. The wound will close, and the AC will remain deep due to the AC maintainer.
- Ensure that the graft is centrally placed over the recipient cornea with the forceps. Inject a small amount of air (2mm size) with an air cannula through the temporal wound, beneath the graft – the air bubble ensures that the graft continues to float up against the recipient stromal surface. Do not attempt to insert a large air bubble at this stage. If the graft is not completely central, a large bubble will quickly press the graft up against the recipient stromal surface, initiating sticking and make subsequent graft positioning difficult.
- Release the graft from the forceps and retract the forceps out of the eye. If there is concern that the wound is not fully self-sealing, and there is high vitreous pressure (in which case the AC may shallow and the graft may shoot out through the wound if it is released) then before releasing the graft, the AC maintainer may be shut off or reduced, and the lip of the wound may be held down with forceps before releasing the graft fully. Full wound closure, air tamponade and completion of the DMEK procedure should then be performed in the usual manner.

**Sterilisation**

- The EndoGlide is a **SINGLE USE ONLY** device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO)
- DO NOT RE-STERILISE.** Network Medical Products Ltd. assumes no liability for devices that have been re-sterilised by healthcare facilities.

**Hazards associated with RE-USE of this Single Use Only device:**

- This single use device has not been validated for re-use. If you re-use a device you may be held **Legally Liable** for the safe performance.
- Cross-contamination and infection risks to patients, including the transmission of:
  - CJD & Variant CJD.
  - Prion Diseases.
  - Bacterial Endotoxins.
  - Hepatitis B & Hepatitis C.
  - Risks posed by HIV and AIDS

Note : Reprocessing may not remove all viable micro-organisms. Abnormal proteins associated with prion diseases eg. Creutzfeldt-Jakob disease (CJD) and variant Creutzfeldt-Jakob disease, have been identified in corneal tissue and are very resistant to all conventional methods of decontamination.
- Device failure through material fatigue or degradation caused by initial use and design. Plastics can be weakened, warped or become brittle. The accuracy, function and performance of the device **will be seriously affected if the device is re-sterilised.**
- Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.
- Separating components** - Forcibly separating the introducer from the cartridge after primary use will lead to leakage of aqueous from the back of the cartridge and loss of intraocular pressure in any re-use.

**Also Recommended for use with the DMEK Endoglide**

|        |                                 |
|--------|---------------------------------|
| 53-951 | Tan EndoGlide Placement Forceps |
| 53-952 | Tan EndoGlide Loading Forceps   |

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**PRODUCT NAME DMEK ENDOGLIDE**

**PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

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