
'PRODUCT NAME: EAR PACK

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GB - INSTRUCTIONS FOR USE

Description

The products in the Network range of PVA sponge products have been designed to meet the requirements of the surgeon where there is a requirement to staunch blood loss after invasive surgery or traumatic injury. The products are quick and easy to use and reduce the time required to create haemostasis where blood loss is critical. The products are available in a variety of designs, widths and lengths to best suit the surgeon's preference and needs. The range is further extended to accommodate for the requirements of children and adults of Asian origin.

Intended Use

PVA sponge ENT products are designed to staunch blood loss after invasive surgery or traumatic injury in ENT surgery. The PVA sponge reacts quickly to absorb blood and body fluids. PVA is biocompatible and in some instances can be used for post-operative tissue support and to carry medication to the operative site (where specified in the Instructions for Use). Post Op Ear Packs are used to prevent stenosis of the ear canal following surgery and will not expand lengthways making it suitable to support a graft following Tympanoplasty.

CAUTIONS

- This device is supplied STERILE and ready to use.
- This device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal Law restricts sales and use to or on the order of a Physician.
- The pack must be kept moist at all times by applying drops 2-3 times daily.
- It is recommended that packs remain in place for no longer than 21 days.

Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

INSTRUCTIONS

Insertion

Lubricate leading edge of pack with antibiotic ointment. Gelfoam or antibiotic ointment may be inserted before the packing is positioned. Grasp the ear pack with forceps.

Once in place hydrate with medication saline or sterile water. A single drop of saline will initially expand the pack if a glycerine based medication is used.

Removal

Always saturate NETCELL ear packs before removal to ensure softness and atraumatic removal. Grasp pack with forceps and gently withdraw.

Sterilisation

This device is supplied sterile by Gamma Irradiation and are designated for SINGLE USE ONLY.

'PRODUCT NAME: EAR PACK

'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES

1. Single use devices have not been validated for re-use.
If you re-use a device you may be held Legally Liable for the safe performance.
2. Cross-contamination and infection risks to patients . Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials

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