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**'PRODUCT NAME: SERIES 5000™ NASAL PACK**  
**'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

			<p>Stati Uniti d'America limita la vendita di questo dispositivo ai medici o su prescrizione del medico</p> <p>NL - Waarschuwing: Volgens de Amerikaanse Federale wetgeving is verkoop van dit apparaat uitsluitend toegestaan door of op voorschrift van een arts.</p> <p>NO - Forsiktig: I henhold til føderale lover i USA skal dette utstyret bare selges av lege eller etter ordinasjon fra lege</p> <p>PT – Cuidado: A Legislação Federal dos EUA restringe a venda deste dispositivo a médicos ou mediante prescrição médica</p> <p>SV - Försiktighet: I USA får enligt federal lag denna produkt endast försäljas av läkare eller enligt läkares ordination</p>				
--	--	--	---	--	--	--	--

**GB - INSTRUCTIONS FOR USE**

**Description**

The products in the Network range of PVA sponge products have been designed to meet the requirements of the surgeon where there is a requirement to staunch blood loss after invasive surgery or traumatic injury. The products are quick and easy to use and reduce the time required to create haemostasis where blood loss is critical. The products are available in a variety of designs, widths and lengths to best suit the surgeon's preference and needs. The range is further extended to accommodate for the requirements of children and adults of Asian origin.

**Intended Use**

PVA sponge ENT products are designed to staunch blood loss after invasive surgery or traumatic injury in ENT surgery. The PVA sponge reacts quickly to absorb blood and body fluids. PVA is biocompatible and in some instances can be used for post-operative tissue support and to carry medication to the operative site (where specified in the Instructions for Use). Series 5000 Nasal packs are single hydration packs that will remain moist for up to 7 days and requires no hydration prior to removal. The external pouch provides non-stick interface against tissue ingrowth which facilitates much easier removal and patient comfort.

**CAUTIONS**

- This device is supplied STERILE and ready to use.
- This device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
- Do not leave packing in for more than 4 days, or serious side effects may occur.
- Always saturate packing before removal.
- String to be used as an indicator only.
- Securing the string to the patient's cheek will prevent the pack moving backwards, thus eliminating the possibility of the pack being imbibed further and causing an obstruction.

**POSSIBLE ADVERSE EFFECTS**

Tampons are associated with toxic shock syndrome (TSS). In case of symptoms such as sudden fever, vomiting, diarrhoea, dizziness, fainting or near fainting and for a rash appearing like a sunburn the treatment with NETCELL has to be interrupted immediately.

**Incident Reporting**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

**INSTRUCTIONS**

**Insertion**

Apply antibiotic ointment to all surfaces of the pack.

Hold the packing with forceps or fingers and insert gently into the desired packing area. The sponge portion of the pack should be within the nose and the neck portion of the pouch should remain exposed.

Once in place tilt the patient's head back. Open the plastic neck portion of the pouch between thumb and forefinger, insert the tip of a syringe and gradually irrigate with a minimum of 10cc's of saline to ensure full expansion.

Secure the pack in place by taping the string to the patient's cheek. If both sides of the nose are packed, the strings may be tied together.

**Removal**

Grasp the packing and gently remove.

Removal to be carried out by a physician.

**Sterilisation**

This device is supplied sterile by Gamma Irradiation and are designated for SINGLE USE ONLY.

**HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES**

1. Single use devices have not been validated for re-use.
  - If you re-use a device you may be held Legally Liable for the safe performance.
2. Cross-contamination and infection risks to patients . Including transmission of:
  - CJD & Variant CJD.
  - Prion Diseases.
  - Bacterial Endotoxins.
  - Hepatitis B & Hepatitis C.
  - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials

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**'PRODUCT NAME: SERIES 5000™ NASAL PACK**  
**'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

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**Network Medical Products Ltd**
Advena Limited