


PRODUCT NAME PVA AND CELLULOSE PRODUCTS FOR OPHTHALMIC APPLICATIONS
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

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GB - INSTRUCTIONS FOR USE

Description

The NETWORK MEDICAL range is a complete line of sponge products used for the management of fluids during ophthalmic procedures. Our PVA and Cellulose sponge products are used during cataract and other refractive surgeries.

The products in the Network range of PVA and Cellulose sponge products have been designed to meet the requirements of the surgeon where there is a requirement to staunch blood loss after invasive surgery or traumatic injury. The products are quick and easy to use and reduce the time required to create haemostasis where blood loss is critical.

Intended Use

PVA and Cellulose sponge ophthalmic products are designed to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. The PVA and Cellulose sponge reacts quickly to absorb blood and body fluids. PVA and Cellulose sponge is biocompatible.

PVA instrument wipes are used to wipe down instruments after use.

The headbands are made from PVA and are used to absorb perspiration during procedures. They are worn by the operating staff and not the patient.

CAUTIONS

- This device is supplied STERILE and ready to use.
- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal Law restricts sales and use to or on the order of a Physician

Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

'PRODUCT NAME PVA AND CELLULOSE PRODUCTS FOR OPHTHALMIC APPLICATIONS
'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

Sterilisation

The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Irradiation.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES

- Single use devices have not been validated for re-use.
If you re-use a device you may be held Legally Liable for the safe performance.
- Cross-contamination and infection risks to patients . Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
- Device failure through material fatigue or degradation caused by initial use and design.
- Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials

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