PRODUCT NAME PVA AND CELLULOSE PRODUCTS FOR OPHTHALMIC APPLICATIONS
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

Symbols Used

<p>LOT</p> <p>GB - Lot Number CS - Číslo šarže DE - Chargennummer ES - Número de lote FR - Numéro de lot IT - Numero di lotto PT - Número de lote SV - Partnummer</p>	<p></p> <p>GB - Use Until Date CS - Spotřebujte do DE - Verwendbar bis Datum ES - Usar antes de la fecha FR - Date limite d'utilisation IT - Utilizzare entro PT - Utilizar até SV - Används före</p>	<p></p> <p>GB - Do Not Re-use CS - Nepoužívejte opakovaně DE - Nicht wiederverwenden ES - No reutilizar FR - Ne pas réutiliser IT - Non riutilizzare PT - Não reutilizar SV - Får ej återanvändas</p>	<p></p> <p>GB - Do not use if product is opened or damaged. CS - Pokud je výrobek otevřený nebo poškozený, nepoužívejte ho. DE - Nicht verwenden, wenn das Produkt bereits offen oder beschädigt ist. ES - No utilizar si el producto está abierto o dañado. FR - Ne pas utiliser si l'emballage est ouvert ou endommagé IT - Non usare se la confezione è aperta o danneggiata PT - Não utilizar se o produto estiver aberto ou danificado SV - Använd inte om produkten har öppnats eller skadats.</p>	<p></p> <p>GB - Do not re-sterilise CS - Nesterilizujte opakovaně DE - Nicht erneut sterilisieren ES - No reesterilizar FR - Ne pas restériliser IT - Non risterilizzare PT - Não voltar a esterilizar SV - Får ej återsteriliseras</p>	<p>STERILE R</p> <p>GB - Sterilised by Irradiation CS - Sterilizováno ozářením DE - sterilisiert durch Bestrahlung ES - Esterilizado por irradiación FR - Stérilisé par irradiation IT - Sterilizzato mediante irradiazione PT - esterilizado por irradiação SV - Steriliserad med strålning</p> <p></p> <p>GB - Non Sterile CS - Nesterilní DE - Nicht-steril ES - No estéril FR - Non stérile IT - Non Sterile PT - Não estéril SV - Icke-steril</p>	<p></p> <p>GB - Manufacturer CS - Výrobce DE - Hersteller ES - Fabricante FR - Fabricant IT - Fabbicante PT - Fabricante SV - Tillverkare</p>	<p>EC REP</p> <p>GB - EU Authorised Representative CS - Autorizovaný zástupce v EU DE - EU-Bevollmächtigter ES - Representante autorizado de la UE FR - Représentant UE autorisé IT - Mandatario nell'Unione europea PT - Representante autorizado na UE SV - Auktoriserad representant i EU.</p>	<p></p> <p>GB - Date of Manufacture CS - Datum výroby DE - Herstellungsdatum ES - Fecha de fabricación FR - Date de fabrication IT - Data di fabbricazione PT - Data de fabrico SV - Tillverkningsdatum</p>	<p></p> <p>GB - Keep Dry CS - Udržujte v suchu DE - Trocken aufbewahren ES - Mantener seco FR - Maintenir au sec IT - Conservare a secco PT - Manter seco SV - Förvaras torr</p>
<p>REF</p> <p>GB - Catalogue Number CS - Katalogové číslo DE - Bestellnummer ES - Número de catálogo FR - Numéro du catalogue IT - Riferimento di catalogo PT - Número de catálogo SV - Katalognummer</p>	<p></p> <p>GB - Caution CS - Upozornění DE - Vorsicht ES - Advertencia FR - Avertissement IT - Attenzione PT - Advertência SV - Försiktighet</p>	<p>NEC_009</p> <p></p> <p>GB - Consult Electronic Instructions for Use CS - Viz elektronický návod k použití DE - Elektronische Gebrauchsanweisung beachten ES - Consultar instrucciones de empleo electrónicas FR - Consulter le mode d'emploi électronique IT - Consultare le Istruzioni per l'uso elettroniche PT - Consultar as instruções de utilização eletrónicas SV - Se bruksanvisningen i elektroniskt format.</p>	<p>Rx Only</p> <p>GB - Caution: US Federal law restricts this device to sale by or on the order of a physician CS - Upozornění: Federální zákon USA omezuje prodej tohoto prostředku na lékaře nebo jejich objednávku DE - Vorsicht: Nach dem US-amerikanischen Bundesgesetz darf dieses Produkt nur von medizinischen Fachkreisen oder auf Anordnung dieser gekauft werden. ES - Advertencia: la legislación federal de EE. UU. establece la restricción de que este dispositivo debe venderse solo a médicos o por orden de estos FR - Avertissement : La loi fédérale américaine limite la vente et l'utilisation aux médecins ou à la demande d'un médecin IT - Attenzione: la legge federale degli Stati Uniti d'America limita la vendita di questo dispositivo ai medici o su prescrizione del medico. PT - Cuidado: A Legislação Federal dos EUA restringe a venda deste dispositivo a médicos ou mediante prescrição médica. SV - Försiktighet: I USA får enligt federal lag denna produkt endast försäljas av läkare eller enligt läkares ordination.</p>	<p>MD</p> <p>GB - Medical Device CS - Zdravotnický prostředek DE - Medizinprodukt ES - Dispositivo médico FR - Dispositif médical IT - Dispositivo medico PT - Dispositivo médico SV - Medicinteknisk produkt</p>	<p>UDI</p> <p>GB - Unique Device Identifier CS - Jedinečný identifikátor prostředku DE - Eindeutige Medizinprodukt-t-Kennung ES - Identificador único del dispositivo FR - Identifiant unique du dispositif IT - Identificatore univoco del dispositivo PT - Identificador de dispositivo único SV - Unik produktidentifera</p>	<p></p> <p>GB - Sterile barrier system CS - Systém sterilní bariéry DE - Steriles Barriersystem ES - Sistema de barrera estéril FR - Système de barrière stérile IT - Sistema di barriera sterile PT - Sistema de barreira estéril SV - Sterilbarriarsystem</p>	<p></p> <p>GB - Sterile barrier system with an additional packaging layer inside CS - Systém sterilní bariéry s další vrstvou obalu uvnitř DE - Steriles Barriersystem mit zusätzlicher Schicht in der Verpackung ES - Sistema de barrera estéril con una capa de embalaje adicional en el interior FR - Système de barrière stérile avec emballage protecteur supplémentaire à l'intérieur IT - Sistema di barriera sterile con uno strato di confezionamento aggiuntivo al suo interno PT - Sistema de barreira estéril com uma camada interior de embalagem adicional. SV - Sterilbarriarsystem med ett extra förpackningslager inuti</p>	<p></p> <p>GB - Non-sterile packaging with sterile barrier inside CS - Nesterilní balení se sterilní bariérou uvnitř DE - Nicht-sterile Verpackung mit Sterilbarriersystem ES - Envase no estéril con barrera estéril en el interior FR - Emballage non stérile contenant une protection stérile IT - Confezione non sterile con barriera sterile interna PT - Embalagem não estéril com barreira estéril no interior SV - Icke-steril emballage med sterilbarriär inuti</p>	<p></p> <p>GB - Keep away from sunlight CS - Chraňte před slunečním zářením DE - Vor Lichteinstrahlung schützen ES - Mantener alejado de la luz solar FR - Conservé à l'abri de la lumière IT - Tenere al riparo della luce solare PT - Manter afastado da luz solar SV - Förvaras skyddat från solljus</p>

PRODUCT NAME PVA AND CELLULOSE PRODUCTS FOR OPHTHALMIC APPLICATIONS
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

This Instructions for use relates to the following products:

GB - This Instructions for use relates to the following products:

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GB - INSTRUCTIONS FOR USE
Product Description:

The PVA and Cellulose Spears, Points, Drains and Wicks are a complete line of sponge products used for the management of fluids during ophthalmic procedures. Our PVA and Cellulose sponge products are used during cataract and other refractive surgeries. The PVA products are manufactured from 100% biocompatible PVA sponge and supplied in various shapes and sizes according to use. The Cellulose sponge products are manufactured from biocompatible Cellulose (approx. 78%) and Pulp (approx. 22%) and supplied in various shapes and sizes according to use.

The PVA Shields are placed on the cornea when it is necessary to moisten the cornea and shield the retina from intense operating light during ophthalmic surgery. These products are manufactured from 100% PVA Sponge in various sizes. Two variants are offered – for general ophthalmic procedures (Eye Shield) and for Lasik Procedures (Lasik Shield).

PVA Instrument Wipes are manufactured from 100% biocompatible PVA sponge in sheet form in various sizes. They are used to help clean medical devices after use.

Intended Population / Patient Target Groups:

Instrument wipes are not intended to be used on patients. All other devices are intended to be used on patients undergoing ophthalmic procedures. Cellulose sponge products however are contraindicated on patients undergoing (LASIK) procedures. There are no other restrictions on patient population. Medical assessment is necessary and due care is to be exercised in cases of eye infections and injuries to the cornea.

Intended Users:

Healthcare professionals only; trained ophthalmic surgeons and ophthalmic theatre staff including trained decontamination professionals.

PRODUCT NAME PVA AND CELLULOSE PRODUCTS FOR OPHTHALMIC APPLICATIONS
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

Intended Purpose:

- **Statement for Intended purpose for PVA spears:**
A body and/or body orifice contact, single-use triangular piece of medical grade, highly absorbent PVA supplied with attachment to a polypropylene handle. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. Can also be used for tissue manipulation. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA points:**
A body and/or body orifice contact, single-use triangular piece of medical grade, highly absorbent PVA. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. Can also be used for tissue manipulation. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for Cellulose spears:**
A body and/or body orifice contact, single-use triangular piece of medical grade, highly absorbent cellulose supplied with attachment to a polypropylene handle. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. Can also be used for tissue manipulation. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. They are contraindicated for patients undergoing LASIK surgery.
- **Statement for Intended purpose for Cellulose points:**
A body and/or body orifice contact, single-use triangular piece of medical grade, highly absorbent cellulose. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. Can also be used for tissue manipulation. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. They are contraindicated for patients undergoing LASIK surgery.
- **Statement for Intended purpose for Cellulose Strips:**
A body or body orifice contact, single-use rectangular strip of medical grade, highly absorbent cellulose used during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. Can also be used for tissue manipulation. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. They are contraindicated for patients undergoing LASIK surgery.
- **Statement for Intended purpose for PVA eye drains:**
A body and/or body orifice contact, single-use rectangular strip of medical grade, highly absorbent PVA which drains into an attached polythene bag which has a self-adhesive strip to permit it to be secured to the drape during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. The PVA portion of the drain has a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA LASIK eye drains:**
A body and/or body orifice contact, single-use piece of medical grade, highly absorbent PVA used during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. They can be used to separate the flap from the cornea along the hinge during LASIK procedures. The PVA LASIK eye drain has a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA wicks:**
A body or body orifice contact, single-use rectangular strip of medical grade, highly absorbent PVA used during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA pledget:**
A body or body orifice contact, single-use rectangular strip of medical grade, highly absorbent PVA used during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA instrument wipes:**
A non-contact, single-use sheet of PVA intended to be used in a sterile condition by healthcare professionals to help clean medical devices prior to full decontamination. The instrument wipe is only to be handled with medical grade gloves and as such is non-contact with a cumulative transient use of less than 60 minutes. There are no known contraindications for this device. These devices are not intended to be used on a patient.
- **Statement for Intended purpose for PVA headbands:**
A body-contact, single-use headband formed from a piece of medical grade, highly absorbent PVA with an elastic strap, intended to be used by healthcare professionals for cooling down their foreheads by hydrating with water and/or absorbing perspiration. The PVA headband is positioned before applying the medical grade gloves and has a cumulative transient use of less than 60 minutes. There are no known contraindications for this device. These devices are not intended to be used on a patient.
- **Statement for Intended purpose for PVA LASIK eye shields:**
A body orifice contact, single-use strip of medical grade and highly absorbent PVA. It is intended to be used by healthcare professionals in a sterile condition on any patient for keeping the surface of the cornea moist and cool during LASIK procedures, and for shielding the retina from intense operating light. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye) and/or the cornea. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA corneal light shields:**
A body orifice contact, single-use circular shaped piece of medical grade and highly absorbent PVA. It is intended to be used by healthcare professionals in a sterile condition on any patient for keeping the surface of the cornea moist and cool during ophthalmic procedures, and for shielding the retina from intense operating light. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye) and/or the cornea. There are no known contraindications for these devices.

Indications for use:

- **PVA spears, PVA points, cellulose spears, cellulose points and cellulose strips:**
Indicated for all microsurgical ophthalmic procedures for tissue manipulation and management of fluids in the operative field. The use of these products is based on their ability to hydrate quickly and efficiently. Used to absorb any blood seeping from the operative site or excess saline used during the procedure. The sponge comes in contact with the tissue but is not invasive of the wound. The design of the sponge enables the surgeon to precisely place the very tip onto the operative site and for hydration to take place quickly from that point. No medications or other substances should be delivered by these products. The cellulose medical devices are not recommended for LASIK procedures.
- **PVA drains and PVA wicks:**
Indicated for all microsurgical ophthalmic procedures for tissue manipulation and management of fluids in the operative field. The wick is placed between the eyeball and the lower eyelid to drain away excess liquid from the operative site or injury. Saline can be used to moisten the wick prior to placement to increase the speed of the wicking.
- **PVA headbands:**
Are used during a procedure in which the user's temperature produces perspiration, not used on the patient.
- **PVA LASIK eye shields:**
Indicated for placement on the cornea when it is necessary to moisten the cornea and shield the retina from intense operating light during ophthalmic surgery.
- **PVA corneal light shields:**
Indicated for placement on the cornea to moisten the cornea and shield the retina from intense operating light during ophthalmic surgery. Usually placed on the cornea during a break in the operating procedure to protect the delicate internal membranes of the eye from damage associated with the heat generated from the microscope light beam. Also maintains moisture on the surface of the cornea.
- **Wicks and wipes:**
A combination product of an instrument wipe and an eye wick. The wick can be detached and placed between the eyeball and the lower eyelid to drain away excess liquid from the operative site or injury. Saline can be used to moisten the wick prior to placement to increase the speed of the wicking. Instrument wipes are used to help clean medical devices.
- **Pledgets:**
Indicated for use when it is necessary to have a continuous collection and control of excess fluids during ophthalmic surgery.

Contraindications:

Cellulose sponge products are contraindicated on patients undergoing LASIK procedures. PVA sponge is a very safe fully biocompatible material and products manufactured from this material are safe for transient and short-term use. There have been no reported incidents of adverse reactions to the material when used correctly and there are no specific clinical contraindications for use of any of the products in general procedures. Instrument wipes are contraindicated for patient use.

PRODUCT NAME PVA AND CELLULOSE PRODUCTS FOR OPHTHALMIC APPLICATIONS
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

Side effects:

Possible side effects include:

- Scratching of the ocular mucous membrane by direct excessive forces from a dry device being applied to the eye
- Irritation if any loose particles are not rinsed away from the treatment site
- Cross infection / contamination from incorrect handling
- Foreign body granuloma if any loose particles are not rinsed away from the treatment site
- Bleb inflammation if any loose particles contaminate in or around a bleb

Clinical benefits, see also Intended Purpose:

Instrument wipes are not used in a clinical application.

- **All patient contact PVA products**
 - Manages fluids and staunches blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery.
 - Can be used for tissue manipulation.
 - Ultra-fast wicking action.
 - Highly absorbent.
 - Supplied compressed but expands immediately upon contact with body fluids, saline or balanced salt solution, wicking 10-15 times its own weight away from the operative site.
 - Supplied sterile, single use only, declared 5-year shelf life.
- **Patient contact PVA spears and PVA points and PVA pledgets:**
 - Suitable for Laser-Assisted In-Situ Keratomileusis (LASIK) procedures.
- **Patient contact PVA eye drains:**
 - Available with either 80cc or 400cc capacity bag which has a self-adhesive strip to permit it to be secured to a drape.
- **Patient contact PVA LASIK eye drains**
 - Can be used to separate the flap from the cornea along the hinge.
- **Patient contact PVA LASIK eye shields:**
 - Supplied sterile, single use only, declared 5-year shelf life.
 - Used to keep the cornea moist and cool during LASIK procedures.
 - Shields retina from intense operating light.
- **Patient contact PVA corneal light shields:**
 - Supplied sterile, single use only, declared 5-year shelf life.
 - Used to keep the cornea moist and cool during ophthalmic procedures.
 - Shields retina from intense operating light.
- **Non-patient contact PVA headbands:**
 - Provides a cooling effect via evaporation once hydrated with water.
 - Reduces contamination to the operating site by controlling and limiting the perspiration from operating staff.
 - Soft and flexible.
 - Supplied sterile, single use only, declared 5-year shelf life.
- **All cellulose spears, cellulose points, cellulose strips:**
 - Manages fluids and staunches blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery.
 - Can be used for tissue manipulation.
 - Constructed from medical grade biocompatible cellulose.
 - Highly absorbent.
 - Supplied sterile, single use only, declared 5-year shelf life.

CAUTIONS:

- These are delicate medical devices and should be handled with care. Do not apply excessive forces to the headband.
- Medical grade (e.g. CE marked) surgical gloves should be used to handle these devices, cross-contamination and infection risks due to unwanted exposure of bacteria, prion, fungus or virus.
- If using sponge holders or forceps to manipulate the devices, ensure these are CE marked for use with PVA and cellulose sponge products.
- PVA and Cellulose Sponge products are to be stored out of direct sunlight and at normal room temperature, or products may expand compromising the pouch seal.
- Do NOT re-sterilise or re-use. Devices are supplied sterile by irradiation ready for use.
- Sterility and performance are not guaranteed on devices that have passed their Use By date or if the packaging or seals have been damaged. If label including use by date is illegible, do not use.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- As with all PVA and cellulose products loose particles may occur, these should be rinsed from the treatment site.
- Do not cut the device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination. These devices should be disposed of after use by following your facilities approved procedures for contaminated waste.
- Control use of devices to ensure none are lost during the procedures, devices do not show up on X-Ray or MRI scanners. Check quantities are correct.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.

Incident Reporting:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Sterilisation:

The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Irradiation.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

- Single use devices have not been validated for re-use. If you re-use a device you may be held Legally Liable for the safe performance.
- Cross-contamination and infection risks due to unwanted exposure of bacteria, prions, fungus or virus.
- Device failure through material fatigue or degradation caused by initial use and design.
- Patient injury from device failure or from residues of decontamination agents absorbed into the materials.

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