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
PRODUCT NAME PVA AND CELLULOSE PRODUCTS FOR OPHTHALMIC APPLICATIONS
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

Lasik Drain 40mm x 14mm x 10.5mm opening	Lasik drén, otvor 40 mm x 14 mm x 10,5 mm	Lasik-Sammelbeutel 40 mm x 14 mm x 10,5 mm Öffnung	Drenaje para cirugía Lasik 40 mm x 14 mm Abertura de 10.5 mm	Drain Lasik 40mm x 14mm x 10, 5mm (diamètre)	Drenaggio Lasik apertura 40 mm x 14 mm x 10, 5 mm	Lasik-drain 40 mm x 14 mm x 10,5 mm opening	Dren Lasik 40 mm x 14 mm x 10,5 mm otwór	Dreno Lasik abertura 40 mm x 14 mm x 10,5 mm	Lasik-drän 40 mm x 14 mm x 10,5 mm öppning
Lasik Drain	Lasik drén	Lasik-Sammelbeutel	Drenaje para cirugía Lasik	Drain Lasik	Drenaggio Lasik	Lasik-drain	Dren Lasik	Dreno Lasik	Lasik-drän
PVA Instrument Wipes	PVA nástrojové utěrky	Instrumentenwischt ücher aus PVA	Toallitas húmedas de PVA para instrumental	Lingettes PVA à instruments	Salviette per strumenti in PVA	PVA- instrumentendo ekjes	Chusteczki z PVA do narzędzi	Toalhetes para instrument os de PVA	PVA- instrumentorkdu kar
Headband	Čelenka	PVA Stirnband	Cinta para la cabeza	Bandeau	Fascia antisudore	Hoofdband	Opaska	Dispositivo para a cabeça	Pannband
Eye Pledget	Oční tampon (obdélníče k)	Augentupfer	Minicompres ocular	Éponge oculaire	Compresa oculare	Oogtampon	Tampony do oczu	Compresa ocular	Ögonkompres

GB - INSTRUCTIONS FOR USE
Product Description:

The PVA and Cellulose Spears, Points, Drains and Wicks are a complete line of sponge products used for the management of fluids during ophthalmic procedures. Our PVA and Cellulose sponge products are used during cataract and other refractive surgeries. The PVA products are manufactured from 100% biocompatible PVA sponge and supplied in various shapes and sizes according to use. The Cellulose sponge products are manufactured from biocompatible Cellulose (approx. 78%) and Pulp (approx. 22%) and supplied in various shapes and sizes according to use.

The PVA Shields are placed on the cornea when it is necessary to moisten the cornea and shield the retina from intense operating light during ophthalmic surgery. These products are manufactured from 100% PVA Sponge in various sizes. Two variants are offered – for general ophthalmic procedures (Eye Shield) and for Lasik Procedures (Lasik Shield).

PVA Instrument Wipes are manufactured from 100% biocompatible PVA sponge in sheet form in various sizes. They are used to help clean medical devices after use.

Intended Population / Patient Target Groups:

Instrument wipes are not intended to be used on patients. All other devices are intended to be used on patients undergoing ophthalmic procedures. Cellulose sponge products however are contraindicated on patients undergoing (LASIK) procedures. There are no other restrictions on patient population. Medical assessment is necessary and due care is to be exercised in cases of eye infections and injuries to the cornea.

Intended Users:

Healthcare professionals only; trained ophthalmic surgeons and ophthalmic theatre staff including trained decontamination professionals.

Intended Purpose:

- **Statement for Intended purpose for PVA spears:**
A body and/or body orifice contact, single-use triangular piece of medical grade, highly absorbent PVA supplied with attachment to a polypropylene handle. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. Can also be used for tissue manipulation. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA points:**
A body and/or body orifice contact, single-use triangular piece of medical grade, highly absorbent PVA. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. Can also be used for tissue manipulation. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for Cellulose spears:**
A body and/or body orifice contact, single-use triangular piece of medical grade, highly absorbent cellulose supplied with attachment to a polypropylene handle. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. Can also be used for tissue manipulation. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. They are contraindicated for patients undergoing LASIK surgery.
- **Statement for Intended purpose for Cellulose points:**
A body and/or body orifice contact, single-use triangular piece of medical grade, highly absorbent cellulose. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. They are contraindicated for patients undergoing LASIK surgery.
- **Statement for Intended purpose for Cellulose Strips:**
A body or body orifice contact, single-use rectangular strip of medical grade, highly absorbent cellulose used during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. Can also be used for tissue manipulation. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. They are contraindicated for patients undergoing LASIK surgery.
- **Statement for Intended purpose for PVA eye drains:**
A body and/or body orifice contact, single-use rectangular strip of medical grade, highly absorbent PVA which drains into an attached polythene bag which has a self-adhesive strip to permit it to be secured to the drape during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. The PVA portion of the drain has a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA LASIK eye drains:**
A body and/or body orifice contact, single-use piece of medical grade, highly absorbent PVA used during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. They can be used to separate the flap from the cornea along the hinge during LASIK procedures. The PVA LASIK eye drain has a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA wicks:**
A body or body orifice contact, single-use rectangular strip of medical grade, highly absorbent PVA used during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA pledget:**
A body or body orifice contact, single-use rectangular strip of medical grade, highly absorbent PVA used during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.
- **Statement for Intended purpose for PVA instrument wipes:**
A non-contact, single-use sheet of PVA intended to be used in a sterile condition by healthcare professionals to help clean medical devices prior to full decontamination. The instrument wipe is only to be handled with medical grade gloves and as such is non-contact with a cumulative transient use of less than 60 minutes. There are no known contraindications for this device. These devices are not intended to be used on a patient.
- **Statement for Intended purpose for PVA headbands:**
A body-contact, single-use headband formed from a piece of medical grade, highly absorbent PVA with an elastic strap, intended to be used by healthcare professionals for cooling down their foreheads by hydrating with water and/or absorbing perspiration. The PVA headband is positioned before applying the medical grade gloves and has a cumulative transient use of less than 60 minutes. There are no known contraindications for this device. These devices are not intended to be used on a patient.
- **Statement for Intended purpose for PVA LASIK eye shields:**

PRODUCT NAME PVA AND CELLULOSE PRODUCTS FOR OPHTHALMIC APPLICATIONS
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

A body orifice contact, single-use strip of medical grade and highly absorbent PVA. It is intended to be used by healthcare professionals in a sterile condition on any patient for keeping the surface of the cornea moist and cool during LASIK procedures, and for shielding the retina from intense operating light. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye) and/or the cornea. There are no known contraindications for these devices.

- **Statement for Intended purpose for PVA corneal light shields:**

A body orifice contact, single-use circular shaped piece of medical grade and highly absorbent PVA. It is intended to be used by healthcare professionals in a sterile condition on any patient for keeping the surface of the cornea moist and cool during ophthalmic procedures, and for shielding the retina from intense operating light. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye) and/or the cornea. There are no known contraindications for these devices.

Indications for use:

- **PVA spears, PVA points, cellulose spears, cellulose points and cellulose strips:**

Indicated for all microsurgical ophthalmic procedures for tissue manipulation and management of fluids in the operative field.

The use of these products is based on their ability to hydrate quickly and efficiently. Used to absorb any blood seeping from the operative site or excess saline used during the procedure. The sponge comes in contact with the tissue but is not invasive of the wound. The design of the sponge enables the surgeon to precisely place the very tip onto the operative site and for hydration to take place quickly from that point. No medications or other substances should be delivered by these products.

The cellulose medical devices are not recommended for LASIK procedures.

- **PVA drains and PVA wicks:**

Indicated for all microsurgical ophthalmic procedures for tissue manipulation and management of fluids in the operative field. The wick is placed between the eyeball and the lower eyelid to drain away excess liquid from the operative site or injury. Saline can be used to moisten the wick prior to placement to increase the speed of the wicking.

- **PVA headbands:**

Are used during a procedure in which the user's temperature produces perspiration, not used on the patient.

- **PVA LASIK eye shields:**

Indicated for placement on the cornea when it is necessary to moisten the cornea and shield the retina from intense operating light during ophthalmic surgery.

- **PVA corneal light shields:**

Indicated for placement on the cornea to moisten the cornea and shield the retina from intense operating light during ophthalmic surgery.

Usually placed on the cornea during a break in the operating procedure to protect the delicate internal membranes of the eye from damage associated with the heat generated from the microscope light beam. Also maintains moisture on the surface of the cornea.

- **Wicks and wipes:**

A combination product of an instrument wipe and an eye wick. The wick can be detached and placed between the eyeball and the lower eyelid to drain away excess liquid from the operative site or injury.

Saline can be used to moisten the wick prior to placement to increase the speed of the wicking. Instrument wipes are used to help clean medical devices.

- **Pledgets:**

Indicated for use when it is necessary to have a continuous collection and control of excess fluids during ophthalmic surgery.

Contraindications:

Cellulose sponge products are contraindicated on patients undergoing LASIK procedures. PVA sponge is a very safe fully biocompatible material and products manufactured from this material are safe for transient and short-term use. There have been no reported incidents of adverse reactions to the material when used correctly and there are no specific clinical contraindications for use of any of the products in general procedures. Instrument wipes are contraindicated for patient use.

Side effects:

Possible side effects include:

- Scratching of the ocular mucous membrane by direct excessive forces from a dry device being applied to the eye
- Irritation if any loose particles are not rinsed away from the treatment site
- Cross infection / contamination from incorrect handling
- Foreign body granuloma if any loose particles are not rinsed away from the treatment site
- Bleb inflammation if any loose particles contaminate in or around a bleb

Clinical benefits, see also Intended Purpose:

Instrument wipes are not used in a clinical application.

- **All patient contact PVA products**
 - Manages fluids and staunches blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery.
 - Can be used for tissue manipulation.
 - Ultra-fast wicking action.
 - Highly absorbent.
 - Supplied compressed but expands immediately upon contact with body fluids, saline or balanced salt solution, wicking 10-15 times its own weight away from the operative site.
 - Supplied sterile, single use only, declared 5-year shelf life.
- **Patient contact PVA eye drains:**
 - Available with either 80cc or 400cc capacity bag which has a self-adhesive strip to permit it to be secured to a drape.
- **Patient contact PVA LASIK eye drains**
 - Can be used to separate the flap from the cornea along the hinge.
- **Patient contact PVA LASIK eye shields:**
 - Supplied sterile, single use only, declared 5-year shelf life.
 - Used to keep the cornea moist and cool during LASIK procedures.
 - Shields retina from intense operating light.
- **Patient contact PVA corneal light shields:**
 - Supplied sterile, single use only, declared 5-year shelf life.
 - Used to keep the cornea moist and cool during ophthalmic procedures.
 - Shields retina from intense operating light.
- **Non-patient contact PVA headbands:**
 - Provides a cooling effect via evaporation once hydrated with water.
 - Reduces contamination to the operating site by controlling and limiting the perspiration from operating staff.
 - Soft and flexible.
 - Supplied sterile, single use only, declared 5-year shelf life.
- **All cellulose spears, cellulose points, cellulose strips:**
 - Manages fluids and staunches blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery.
 - Can be used for tissue manipulation.
 - Constructed from medical grade biocompatible cellulose.
 - Highly absorbent.
 - Supplied sterile, single use only, declared 5-year shelf life.

CAUTIONS:

- These are delicate medical devices and should be handled with care. Do not apply excessive forces to the headband.
- Medical grade (e.g. CE marked) surgical gloves should be used to handle these devices, cross-contamination and infection risks due to unwanted exposure of bacteria, prion, fungus or virus.
- If using sponge holders or forceps to manipulate the devices, ensure these are CE marked for use with PVA and cellulose sponge products.
- PVA and Cellulose Sponge products are to be stored out of direct sunlight and at normal room temperature, or products may expand compromising the pouch seal.
- Do NOT re-sterilise or re-use. Devices are supplied sterile by irradiation ready for use.
- Sterility and performance are not guaranteed on devices that have passed their Use By date or if the packaging or seals have been damaged. If label including use by date is illegible, do not use.

PRODUCT NAME PVA AND CELLULOSE PRODUCTS FOR OPHTHALMIC APPLICATIONS
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- As with all PVA and cellulose products loose particles may occur, these should be rinsed from the treatment site.
- Do not cut the device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination. These devices should be disposed of after use by following your facilities approved procedures for contaminated waste.
- Control use of devices to ensure none are lost during the procedures, devices do not show up on X-Ray or MRI scanners. Check quantities are correct.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.

Incident Reporting:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Sterilisation:

The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Irradiation.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

- Single use devices have not been validated for re-use. If you re-use a device you may be held Legally Liable for the safe performance.
- Cross-contamination and infection risks due to unwanted exposure of bacteria, prions, fungus or virus.
- Device failure through material fatigue or degradation caused by initial use and design.
- Patient injury from device failure or from residues of decontamination agents absorbed into the materials.

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