NETWORK
ENT PRODUCTS

VENTILATION TUBES
AND GROMMETS

PRODUCT INFORMATION
AND INSTRUCTIONS

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NETWORK MEDICAL PRODUCTS LTD

Coronet House, Kearsley Road,
Ripon, North Yorkshire, HG4 2SG, UK

Tel: +44 (0)1765 609555

Fax: +44 (0)1765 608476

e-mail: info@networkmedical.co.uk

EC REP

Advena Limited
Tower Business Centre
2nd Floor, Tower Street, Swatar
BKR 4013, Malta

CE 1639

SYMBOLS

REF	Catalogue Number		Use-by date
	Date of Manufacture	STERILE EO	Sterilized using ethylene oxide
LOT	Batch code		GB - Sterile barrier system with an additional packaging layer inside
	Manufacturer		Do not resterilise
	Do not use if package is damaged and consult instructions for use		Do not re-use
Rx Only	CAUTION: US Federal Law restricts sales and use to or on the order of a Physician		Caution
MD	Medical Device	NEW 010	Consult Electronic instructions for use
EC REP	Authorized representative in the European Community	UDI	Unique Device Identifier

VENTILATION TUBES - INSTRUCTIONS FOR USE

Device Description:

Network Medical ventilation (tympanostomy) tubes are small tubular implants available in a variety of biocompatible materials including silicone elastomer and fluoroplastic. Designs and sizes are available with single or multiple flanges to satisfy various surgical techniques for insertion and to facilitate short or long term communication of the auditory canal with the middle ear.

Functional Characteristics:

The materials used to manufacture the Network Medical ventilation tubes are silicone and PTFE Polymer. These materials have long histories of safe clinical use as implanted medical devices. They are known for their smoothness of texture, inertness, biocompatibility and effectiveness for both short and long-term middle ear drainage and ventilation.

Intended Use:

Ventilation tubes are indicated for the treatment of chronic otitis media with effusion, recurrent acute otitis media, tympanic membrane atelectasis and complications of acute otitis media in children. The tubes which are surgically inserted through an incision made in the tympanic membrane, act as ventilation devices allowing a free exchange of air between the outer ear and the middle ear space equalizing the pressure on both sides of the membrane. When fluid is presenting in the middle ear, the ventilation tube can also act as a drain, allowing fluid to drain from the middle ear space to the external auditory canal.

The shape and design of the tube is important in determining its function. Large inner flanges increase the duration of the middle ear ventilation. Tubes with a large lumen and short length are less likely to occlude but may allow easier passage of water. Short term tubes, with an average duration of ventilation of 10 months are shaped like grommets, are easy to insert and have a low rate of obstruction and permanent perforation. Long term ventilation tubes, with an average duration of ventilation of 30 months, have the opposite characteristics, and are typically T-shaped.

The tubes are therefore available in a variety of designs, sizes and materials to best suit the surgeon's preference and need.

Intended lifetime:

Ventilation tubes and grommets have a maximum expected lifetime of 18 months.

Indications for Use:

The surgeon must use medical judgement and consider the patient's medical history prior to a decision to surgically insert a ventilation tube. Pathologic conditions for which ventilation tubes are indicated include but are not limited to:

- Chronic otitis media with effusion characterised as either serous, mucoid or purulent.
- Recurrent acute otitis media which fails to respond satisfactorily to alternative therapies.
- A patient with a history of high negative middle ear pressure which may be associated with conductive hearing loss, otalgia, vertigo and/or tinnitus.
- Atelectasis resulting from retraction pocket of the tympanic membrane or Eustachian tube dysfunction.

Caution:

- This device is supplied STERILE and ready to use.
- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal Law restricts sales and use to or on the order of a Physician
- Ventilation tubes and grommets supplied with wire contain nickel which may cause allergic reaction in patients with sensitivity to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.

Incident Reporting

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Implant Removal:

Any implant that has to be removed should be removed by gently grasping the rim or wire extrusion, using a pair of micro alligator forceps and pulling gently.

Contraindications:

Some conditions for which ventilation tubes may not be the therapy of choice include

- Cases of otitis media which respond favourably to drug therapy alone
- Cases of otitis media in which, in the medical judgement of the surgeon, a myringotomy alone or tympanocentesis will suffice as treatment for the condition.
- Cases where the patient has severe Eustachian tube dysfunction, Downs Syndrome or Cleft Palate as these patients can have a higher complication rate.

Possible Adverse Effects:

- The middle ear may develop subsequent secondary infections from either water or airborne pathogens. (Care should be taken to avoid water entering the ear).
- Clogging of the tube which may require removal and replacement, and only for use in ears evacuated of effusion.
- Persistent perforation of the tympanum after tube extrusion
- Granulomatous reaction of the tympanum necessitating tube removal
- Premature extrusion of the tube.
- Failure of the tube to self-extrude requiring medical intervention for removal.
- Tympanosclerosis / myringosclerosis
- tympanic membrane atrophy / Membrane erosion
- tympanic membrane retraction and atelectasis
- formation of retraction pockets

Precautions:

- Larger tube diameters and increasing flange size may increase the chances of prolonged tube retention and/or residual perforation.
- It is known that prolonged T-Tube insertion could increase the risk of tympanic membrane perforation. Therefore Network Medical recommends surgical removal of T-Tubes no later than 18 months following surgery.
- In rare circumstances it has been known for tympanostomy tube displacement into the middle ear to occur.

The surgeon is encouraged to carefully choose a tube size, style and material to meet the needs of each specific clinical situation, while minimising adverse effects.

Sterility:

Network Medical ventilation tubes are provided STERILE, being terminally sterilised by Ethylene Oxide (ETO) gas in a highly controlled environment and confirmed through biological testing. The package should be carefully inspected for punctures tears, security of the seal or any other evidence of the sealed pouch having been compromised prior to placement of the contents in the sterile field. This product is sterile in an undamaged, unopened package.

Resterilisation of ventilation tubes is not recommended. Network Medical assumes no liability for products that have been resterilised.

MRI Compatibility:

All products made from Teflon or Silicon are MRI compatible. Any product with steel components may be compromised by the MRI scanner. Please contact Network Medical for further information.

Disposal:

Used devices should be disposed of as clinical waste. If unused, product in original packaging which has exceeded the declared shelf life should be removed from packaging and disposed of as clinical waste.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

1. Single use devices have not been validated for re-use.
If you re-use a device you may be held **Legally Liable** for the safe performance.
2. Cross-contamination and infection risks to patients . Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
 - **Plastics:** Can be weakened, warped or become brittle.
 - **Metals:** Can be damaged or subjected to rusting.
 - **Other materials:** May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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