PRODUCT NAME: PRECISION LASIK BLADE
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

Symbols Used

LOT GB- Lot Number DE - Chargennummer ES - Número de lote FR - Numéro de lot IT - Numero di lotto PT - Número de lote RU — Номер партии	 GB- Use By Date DE – “Verwendbar bis” Datum ES - Usar antes de la fecha FR - Date limite d'utilisation IT - Utilizzare entro PT - Utilizar até RU — Исползовать до	 GB- Do Not Re-use DE - Nicht wiederverwenden ES - No reutilizar FR - Ne pas réutiliser IT - Non riutilizzare PT - Não reutilizar RU — Запрет на повторное применение	 GB- Do not use if product is opened or damaged. DE - Nicht verwenden, wenn das Produkt bereits offen oder beschädigt ist. ES - No utilizar si el producto está abierto o dañado. FR - Ne pas utiliser si l'emballage est ouvert ou endommagé IT - Non usare se la confezione è aperta o danneggiata PT - Não utilizar se o produto estiver aberto ou danificado RU — Не использовать, если изделие вскрыто или повреждено	STERILE EO GB- Sterilised by Ethylene Oxide DE – Ethylenoxid-Sterilisation ES - Estéril por oxido de etileno FR - Stérilisation par oxyde d'éthylène IT - Sterilizzato mediante ossido di etilene PT – Esterilizado por óxido de etileno RU — Стерилизовано этиленоксидом	 GB –Manufacturer DE - Hersteller ES – Fabricante FR - Fabricant IT – Fabbricante PT - Fabricante RU — Изготовитель	 GB – Date of Manufacture DE - Herstellungsdatum ES - Fecha de fabricación FR - Date de fabrication IT - Data di fabbricazione PT - Data de fabrico RU — Дата изготовления	EC REP GB – EU Authorised Representative DE - EU-Bevollmächtigter ES - Representante autorizado de la UE FR - Représentant UE autorisé IT - Mandatario nell'Unione europea PT - Representante autorizado na UE RU — Уполномоченный представитель в Европейском союзе
REF GB- Catalogue Number DE – Bestellnummer ES - Número de catálogo FR- Numéro du catalogue IT - Riferimento di catalogo PT - Número de catálogo RU — Номер по каталогу	 GB- Caution DE – Vorsicht ES – Advertencia FR- Avertissement IT – Attenzione PT - Advertência RU — Осторожно!	NEW_016 GB - Consult Electronic Instructions for Use DE - Elektronische Gebrauchsanweisung beachten FR - Consulter le mode d'emploi électronique ES - Consultar instrucciones de empleo electrónicas IT - Consultare le Istruzioni per l'uso elettroniche PT - Consultar as instruções de utilização eletrónicas RU — Обратитесь к электронной инструкции по применению	R_{only} GB – Caution: US Federal law restricts this device to sale by or on the order of a physician DE – Vorsicht: Nach dem US-amerikanischen Bundesgesetz darf dieses Medizinprodukt nur von medizinischen Fachkreisen oder auf Anordnung dieser gekauft werden ES - Advertencia: la legislación federal de EE. UU. establece la restricción de que este dispositivo debe venderse solo a médicos o por orden de estos. FR - Avertissement : la loi fédérale américaine limite la vente et l'utilisation aux médecins ou à la demande d'un médecin IT - Attenzione: la legge federale degli Stati Uniti d'America limita la vendita di questo dispositivo ai medici o su prescrizione del medico PT – Cuidado: A Legislação Federal dos EUA restringe a venda deste dispositivo a médicos ou mediante prescrição médica. RU — Внимание! В соответствии с федеральным законодательством США изделие может быть продано только врачу или по заказу врача	 GB- Do not re-sterilise DE - Nicht erneut sterilisieren ES - No reesterilizar FR - Ne pas restériliser IT - Non risterilizzare PT - Não voltar a esterilizar RU — Не стерилизовать повторно	MD GB – Medical Device DE – Medizinprodukt ES - Dispositivo medico FR – Dispositif medical IT - Dispositivo medico PT - Dispositivo medico RU — Медицинское изделие	UDI GB – Unique Device Identifier DE - Eindeutige Medizinprodukt-Kennung ES - Identificador único del dispositivo FR – Identifiant unique du dispositif IT - Identificatore univoco del dispositivo PT - Identificador de dispositivo único RU — Уникальный идентификатор устройства	 GB - Sterile barrier system with an additional packaging layer inside DE – Steriles Barriersystem mit zusätzlicher Schicht in der Verpackung ES - Sistema de barrera estéril con una capa de embalaje adicional en el interior FR - Système de barrière stérile avec emballage protecteur supplémentaire à l'intérieur IT - Sistema di barriera sterile con uno strato di confezionamento aggiuntivo al suo interno PT – Sistema de barreira estéril com uma camada interior de embalagem adicional. RU — Система стерильного барьера с дополнительным слоем упаковки внутри

GB - INSTRUCTIONS FOR USE

ES3030	The NETWORK blade is a single use blade designed for use with the Hansatome® Microkeratome.
ES5050	The NETWORK blade is a single use blade designed for use with the Zyoptix® XP Microkeratome.
ES7070	The NETWORK blade is a single use blade designed for use with the Moria® M2 Microkeratome.
ES9090	The NETWORK blade is a single use blade designed for use with the Nidek® MK-200 Microkeratome.

Description

Network Medical Precision LASIK Blades are generic microkeratome blades designed specifically to fit the Bausch & Lomb Hansatome, Zyoptix, Moria and Nidek Microkeratome systems used to create the corneal flap prior to the laser reshaping.

Intended Use

Plano Precision LASIK Blade: The Plano model is used to create a flap that approximates the average flap thickness produced by the given microkeratome head.

Minus Precision LASIK Blade:

A minus model can be used to create a thin flap and maximise the amount of residual stroma for cases of higher correction. Patients with thicker corneas often experience thick flaps because of the cornea compression that occurs during the microkeratome pass. To adjust for the normally thicker flap on these patients, a Minus blade can be used to compensate for the thicker stroma to obtain a desired flap thickness.

Plus Precision LASIK Blade: The Plus model can be used when there is a need to create a thicker flap than the microkeratome head usually creates. This can be especially beneficial in cases to cut below an original flap depth to re-treat a patient.

CAUTIONS:

- Visually inspect blade for handling damage prior to installing in microkeratome.
- Following the Microkeratome Equipment Manual, carefully load the blade into the Microkeratome head. The blade should load easily, without excessive force. After loading the blade into the microkeratome, confirm that the blade and holder move smoothly back and forth within the cavity before using. If it does not load easily or move free of resistance, remove the blade and try another.
- The blade is designed to promote operator safety when handled correctly.
- This device is supplied STERILE and ready to use.
- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician

PRODUCT NAME: PRECISION LASIK BLADE

PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

- This product contains nickel which may cause allergic reaction in patients with sensitivity to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
- The device is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface, check the blade integrity under the microscope before use and discard if the blade is damaged.
- For safe disposal, the device should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.

Possible Adverse Effects

Complications which may occur include, but are not limited to keraectasia, punctal keratitis, macular haemorrhages, macular holes, rhegmatogenous retinal detachment, retinal holes, erosion, flap melting/dislocation. Patients with monovision may experience a reduction in binocular vision.

Incident Reporting

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Sterilisation

The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

1. Single use devices have not been validated for re-use.
If you re-use a device you may be held **Legally Liable** for the safe performance.
2. Cross-contamination and infection risks to patients. Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
 - **Plastics:** Can be weakened, warped or become brittle.
 - **Metals:** Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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