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
PRODUCT NAME: PRECISION LASIK BLADE
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

GB - INSTRUCTIONS FOR USE

ES3030	The NETWORK blade is a single use blade designed for use with the Hansatome* Microkeratome.
ES5050	The NETWORK blade is a single use blade designed for use with the Zyoptix* XP Microkeratome.
ES7070	The NETWORK blade is a single use blade designed for use with the Moria* M2 Microkeratome.
ES9090	The NETWORK blade is a single use blade designed for use with the Nidek* MK-200 Microkeratome.

Description

Network Medical Precision LASIK Blades are generic microkeratome blades designed specifically to fit the Bausch & Lomb Hansatome, Zyoptix, Moria and Nidek Microkeratome systems used to create the corneal flap prior to the laser reshaping.

Intended Use

Plano Precision LASIK Blade: The Plano model is used to create a flap that approximates the average flap thickness produced by the given microkeratome head.

Minus Precision LASIK Blade:

A minus model can be used to create a thin flap and maximise the amount of residual stroma for cases of higher correction. Patients with thicker corneas often experience thick flaps because of the cornea compression that occurs during the microkeratome pass. To adjust for the normally thicker flap on these patients, a Minus blade can be used to compensate for the thicker stroma to obtain a desired flap thickness.

Plus Precision LASIK Blade: The Plus model can be used when there is a need to create a thicker flap than the microkeratome head usually creates. This can be especially beneficial in cases to cut below an original flap depth to re-treat a patient.

CAUTIONS:

- Visually inspect blade for handling damage prior to installing in microkeratome.
- Following the Microkeratome Equipment Manual, carefully load the blade into the Microkeratome head. The blade should load easily, without excessive force. After loading the blade into the microkeratome, confirm that the blade and holder move smoothly back and forth within the cavity before using. If it does not load easily or move free of resistance, remove the blade and try another.
- The blade is designed to promote operator safety when handled correctly.
- This device is supplied STERILE and ready to use.
- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician
- The device is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface, check the blade integrity under the microscope before use and discard if the blade is damaged.
- For safe disposal, the device should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.

Possible Adverse Effects

Complications which may occur include, but are not limited to keraectasia, punctal keratitis, macular haemorrhages, macular holes, rhegmatogenous retinal detachment, retinal holes, erosion, flap melting/dislocation. Patients with monovision may experience a reduction in binocular vision.

Incident Reporting

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Sterilisation

The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

1. Single use devices have not been validated for re-use.
If you re-use a device you may be held **Legally Liable** for the safe performance.
2. Cross-contamination and infection risks to patients. Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
 - **Plastics:** Can be weakened, warped or become brittle.
 - **Metals:** Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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PRODUCT NAME: PRECISION LASIK BLADE

PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

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