**PRODUCT NAME: INTERNAL NASAL SPLINTS**  
**PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

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**GB - INSTRUCTIONS FOR USE**

**Description**

Network Silicone Nasal Splints (with integral airway) and Internal Bivalve Splints provide septal support and allow nasal breathing post-op through integral airway.

**Intended Use**

Silicone nasal splints with an integral airway product have a pre-punched hole in anterior tip for easier suturing to prevent dislodgement and aid in septal stabilization.

NETWORK® Bivalve Splints are split for atraumatic insertion and removal and are supplied sterile ready for use. 100% fluoroplastic to prevent adhesions between the septum and lateral wall. They are easy to insert and remove and are pre-punched with suture holes.

**CAUTIONS**

- This device is supplied STERILE and ready to use.
- This device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.

**Incident Reporting**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

**Sterilisation**

The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).

**HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES**

- Single use devices have not been validated for re-use. If you re-use a device you may be held Legally Liable for the safe performance.
- Cross-contamination and infection risks to patients. Including transmission of:
  - CJD & Variant CJD.
  - Prion Diseases.
  - Bacterial Endotoxins.
  - Hepatitis B & Hepatitis C.
  - Risks posed by HIV and AIDS
- Device failure through material fatigue or degradation caused by initial use and design:
  - Other materials:** May degrade, becoming unacceptable when compared to original manufacturing criteria.
- Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

**'PRODUCT NAME: INTERNAL NASAL SPLINTS**  
**'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

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