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
**GB - INSTRUCTIONS FOR USE**
**Description**

NETWORK ENT® Silicone Airway Splints and Fluoroplastic Bivalve Splints are used to provide temporary, postoperative splinting of the nasal septum. These splints are designed to reduce postoperative complication rates for patients undergoing septal surgeries.

**Intended Use/Indications for Use**

NETWORK ENT® Silicone Airway Splints and Fluoroplastic Bivalve Splints are intended to be used postoperatively on septal surgery patients to reduce complication rates, such as adhesion, hematoma, and synechia formation.

**Intended User**

NETWORK ENT® Silicone Airway Splints and Fluoroplastic Bivalve Splints are intended to be used by healthcare professionals, such as, ENT, plastic and cosmetic surgeons.

**WARNINGS AND PRECAUTIONS**

- This device is supplied STERILE and ready to use.
- This device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- Internal nasal splints are intended to be sutured to the nasal septum to prevent splint dislodgement.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.

**Instructions for Application and Removal**

Insert one internal nasal splint into each of the patient's nasal cavities

Suture the nasal splint to the nasal septal wall

Remove sutures and splint in two to ten days

**Disposal**

Dispose of device as clinical waste in line with local protocol.

**Incident Reporting**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

**Sterilisation**

The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).

**HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES**

1. Single use devices have not been validated for re-use. If you re-use a device you may be held Legally Liable for the safe performance.
2. Cross-contamination and infection risks to patients. Including transmission of:
  - CJD & Variant CJD.
  - Prion Diseases.
  - Bacterial Endotoxins.
  - Hepatitis B & Hepatitis C.
  - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
  - **Other materials:** May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

**PRODUCT NAME: INTERNAL NASAL SPLINTS**  
**PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

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