



NETWork

Medical Products



ENT Products



PVA Products



NETCELL NASAL PACK

- For use after septal, turbinate and rhinoplastic surgery.
- Provides gentle, equal pressure within the nasal cavity when absorbing post-operative bleeding.
- Patient comfort is enhanced with the addition of the optional airway tube.
- Airway tube is positioned along the nasal floor. This allows for improved venting and reduces the risk of obstruction.
- Slimpack developed for use in conjunction with a nasal splint following septal surgery.



NETCELL NASAL PACK TUBED

The Netcell Nasal Packing can be supplied with an integral airway tube. The airway tube benefits the patient through easier breathing and improved patient comfort.

Supplied Sterile

Product	Dimensions L x Th x H (cm)	Pack	Code
3.0 cm Nasal Pack	3.0 x 1.5 x 2.0 cm	Box of 20	10-095
4.5 cm Nasal Pack	4.5 x 1.5 x 2.0 cm	Box of 20	10-100
4.5 cm Nasal Slimpack	4.5 x 1.0 x 2.0 cm	Box of 20	10-101
4.5 cm Nasal Slimpack with Tube	4.5 x 1.0 x 2.0 cm	Box of 20	10-101T
6.0 cm Nasal Pack	6.0 x 1.5 x 2.0 cm	Box of 20	10-105
6.0 cm Nasal Pack with Tube	6.0 x 1.5 x 2.0 cm	Box of 20	10-105T
6.0 cm Nasal Slimpack	6.0 x 1.0 x 2.0 cm	Box of 20	10-106
8.0 cm Nasal Pack	8.0 x 1.5 x 2.0 cm	Box of 20	10-110
8.0 cm Nasal Pack with Tube	8.0 x 1.5 x 2.0 cm	Box of 20	10-110T
8.0 cm Nasal Slimpack	8.0 x 1.0 x 2.0 cm	Box of 20	10-111
8.0 cm Nasal Slimpack with Tube	8.0 x 1.0 x 2.0 cm	Box of 20	10-111T
10.0cm Nasal Pack	10.0 x 1.5 x 2.0 cm	Box of 20	10-115



NETCELL ANATOMICAL PACK

- Clinically designed to fill the nasal cavity - providing gentle and equal support and allowing the efficient absorption of post operative bleeding.
- The optional airway tube will enhance patient comfort.

Supplied Sterile

Product	Dimensions L x Th x H	Pack	Code
Anatomical Pack	8.0 x 1.5 x 3.0 cm	Box of 10	10-120
Anatomical Pack with Tube	8.0 x 1.5 x 3.0 cm	Box of 10	10-120T
Anatomical Slimpack	8.0 x 1.0 x 3.0 cm	Box of 10	10-125
Anatomical Slimpack with Tube	8.0 x 1.0 x 3.0 cm	Box of 10	10-125T
Anatomical Pack	5.5 x 1.5 x 2.5 cm	Box of 10	10-126



NETCELL EPISTAXIS PACK

- Immediate, effective control of bleeding.
- Fluid absorption causes the packing to expand producing a tamponading effect.
- The Tubed Epistaxis pack benefits the patient through balanced air pressure and encourages easier breathing. Suction can be achieved through the tube.
- Long Epistaxis for the simple quick packing of the entire nasal cavity, the most efficient way to contain and restrict blood loss.
- Short Epistaxis Pack for the quick convenient treatment of anterior Epistaxis.
- Layered pack allows the packing to be adjusted to the individual patient.

Supplied Sterile

Product	Dimensions L x Th x H	Pack	Code
Epistaxis Pack Long	10.0 x 1.5 x 2.5 cm	Box of 10	10-140
Epistaxis Long with Tube	10.0 x 1.5 x 2.5 cm	Box of 10	10-140T
Epistaxis Slimpack Long	10.0 x 1.0 x 2.5 cm	Box of 10	10-141
Epistaxis Slimpack Long Layered	10.0 x 1.0 x 2.5 cm	Box of 10	10-141L
Epistaxis Pack Medium	8.0 x 1.5 x 2.5 cm	Box of 10	10-145
Epistaxis Slimpack Medium	8.0 x 1.0 x 2.5 cm	Box of 10	10-146
Epistaxis Pack Short	5.5 x 1.5 x 2.5 cm	Box of 20	10-150
Epistaxis Pack Short Layered	5.5 x 1.5 x 2.5 cm	Box of 20	10-150L
Naso-Pak	2.8 x 0.7 x 2.5 cm	Box of 10	10-160




SERIES 5000 SINGLE HYDRATION NASAL PACK

- **No need to rehydrate before removal.**
- This nasal packing is the most effective, user-friendly nasal packing with the benefits of single hydration.
- Patient pain and discomfort is radically reduced.
- Double compressed for easy insertion, the outer covering will accept water based or oil based medication, to the surgeon's preference.

Supplied Sterile

Product	Dimensions L x Th x H	Pack	Code
Single Hydration Nasal Pack	6.0 x 1.5 x 2.0 cm	Box of 10	15-105
Single Hydration Nasal Pack	8.0 x 1.5 x 2.0 cm	Box of 10	15-110
Single Hydration Nasal Slimpack	8.0 x 1.0 x 2.0 cm	Box of 10	15-111



NETCELL SINUS PACK

- Developed as a post-operative packing for F.E.S.S.
- Designed to prevent lateralization of the middle turbinates post-operatively.
- Can be used to pack the Osteomeatal complex to leave the nasal passage clear for easier breathing.

Supplied Sterile

Product	Dimensions L x Th x H	Pack	Code
2.5 cm Sinus Pack Standard	2.5 x 1.20 x 2.0 cm	Pack of 2/box of 10 Packs	20-200
2.5 cm Sinus Pack Small	2.5 x 0.60 x 2.0 cm	Pack of 2/box of 10 Packs	20-202
3.5 cm Sinus Pack Standard	3.5 x 1.20 x 1.2 cm	Pack of 2/box of 10 Packs	20-210
3.5 cm Sinus Pack Medium	3.5 x 0.90 x 1.2 cm	Pack of 2/box of 10 Packs	20-211
3.5 cm Sinus Pack Small	3.5 x 0.60 x 1.2 cm	Pack of 2/box of 10 Packs	20-212



NETCELL EAR-WICK

- Ideal for delivering antibiotics when treating Otitis Externa.
- Instant patient relief when hydrated with medicated otic solution.

Supplied Sterile

Product	Dimensions	Pack	Code
15 mm Ear wick	9 mm dia x 15 mm	1 per pack / box of 50 packs	30-300
12 mm Ear wick	7 mm dia x 12 mm	1 per pack / box of 10 packs	30-301
15 mm Ear wick Fenestrated	9 mm dia x 15 mm	1 per pack / box of 50 packs	30-305
12 mm Ear wick Fenestrated	7 mm dia x 12 mm	1 per pack / box of 10 packs	30-306
24 mm Ear wick	9 mm dia x 24 mm	1 per pack / box of 10 packs	30-310
20 mm Ear wick	7 mm dia x 20 mm	1 per pack / box of 10 packs	30-311



NETCELL POST-OP EAR PACK

- Designed to prevent stenosis of the canal following surgery.
- Designed not to expand lengthwise and is suitable to support a graft following tympanoplasty.

Supplied Sterile

Product	Dimensions	Pack	Code
Ear Pack (Strung)	9 mm dia x 15 mm	1 per pack / box of 10 packs	30-320
Ear Pack (Strung)	12 mm dia x 15 mm	1 per pack / box of 10 packs	30-330
Ear Pack without string	12 mm dia x 24 mm	1 per pack / box of 10 packs	30-340




NETCELL MICROSPHERES

- Developed to control haemostasis in any micro surgical technique.
- The microspheres are hydrated with the appropriate medication and applied to the areas of bleeding.

Supplied Sterile

Product	Dimensions	Pack	Code
Microsphere Small	3 mm dia	Pack of 10 / box of 10 packs	30-380
Microsphere Medium	4 mm dia	Pack of 10 / box of 10 packs	30-381
Microsphere Large	5 mm dia	Pack of 10 / box of 10 packs	30-382
Microsphere Combi Pak	3mm/5mm	50 each size	30-383



INSTRUMENT WIPES

Micro-surgical instruments are expensive and easy to damage. Our Instrument Wipes have been especially developed for the maintenance and care of these delicate surgical products. The wipes are soft but very strong, they will not snag and will easily remove the protein that weakens and affects the surface and edges of instruments. The life of a surgical instrument will be extended if it is cleaned thoroughly immediately after use. The wipes are supplied sterile for use during surgery.




HEADBANDS

Beneficial for surgeons operating in demanding situations, when moistened with water the headband will have an immediate cooling effect through evaporation.

Supplied Non-Sterile

Product	Sterile / Non Sterile	Pack	Code
Instrument Wipe	Sterile	Pack of 1/box of 20	35-900
Headband	Non Sterile	Pack of 1/box of 10	35-950



DCR & LACHRYMAL INTUBATION

- Supplied sterile for single use only
- Sterilisation by Gamma Irradiation
- Declared 5 year shelf life



RETRIEVAL DEVICE FOR CRAWFORD INTUBATION SET

Supplied Sterile

Probe Length	Overall Length	Pack	Code
77 mm	105 mm	Box of 5	51-917



CRAWFORD-TYPE INTUBATION SET

Supplied Sterile

Probe Length	Silicone Tube Length	Pack	Code
150 mm	320 mm	Single pack	51-918

Supplied complete with retrieval device



CANALICULUS INTUBATION SET

Supplied Sterile

Probe Length	Silicone Tube Length	Pack	Code
185 mm	300 mm	Single pack	51-920



DCR BODKIN STRAIGHT

Supplied Sterile


Bodkin Length	Silicone Rod Length	Pack	Code
45 mm / 19 gauge	400 mm	Box of 3	51-922



DCR BODKIN ANGLED

Supplied Sterile

Bodkin Length	Silicone Rod Length	Pack	Code
45 mm / 19 gauge	400 mm	Box of 3	51-923



STAPES PISTON

Trim-to length Teflon (fluoroplastic) stapes piston

- Fluoroplastic shaft may be trimmed to length easily using a sharp knife.
- Loop can be opened and the plastic memory of the fluoroplastic provides secure closure.
- Supplied sterile in single packs.

Supplied Sterile

Shaft Diameter	Loop Diameter	Length	Pack	Code
0.4 mm	0.6 mm	6.0 mm	Single Pack	71-400
0.4 mm	0.6 mm	5.0 mm	Single Pack	71-401
0.4 mm	0.6 mm	4.5 mm	Single Pack	71-402
0.4 mm	0.6 mm	4.0 mm	Single Pack	71-403
0.4 mm	0.6 mm	3.5 mm	Single Pack	71-404
0.6 mm	0.6 mm	6.0 mm	Single Pack	71-405
0.6 mm	0.6 mm	5.0 mm	Single Pack	71-406
0.6 mm	0.6 mm	4.5 mm	Single Pack	71-407
0.6 mm	0.6 mm	4.0 mm	Single Pack	71-408
0.6 mm	0.6 mm	3.5 mm	Single Pack	71-409
0.8 mm	0.6 mm	6.0 mm	Single Pack	71-410
0.8 mm	0.6 mm	5.0 mm	Single Pack	71-411
0.8 mm	0.6 mm	4.5 mm	Single Pack	71-412
0.8 mm	0.6 mm	4.0 mm	Single Pack	71-413
0.8 mm	0.6 mm	3.5 mm	Single Pack	71-414



EPISTAXIS CATHETER

- Rapid and effective nasal bleeding control.
- Streamlined design assists catheter insertion and removal, improving patient comfort.
- All silicone construction ensures the cuffs will remain soft and compliant.
- Minimised adherence to the nasal mucosa.
- Airway permits nasal breathing and access for suctioning.
- Independent inflatable cuffs facilitate proper placement and will conform to nasal anatomy.
- The most effective control of posterior epistaxis.

Supplied Sterile

Description	Pack	Code
Large Epistaxis Balloon Catheter	Box of 1	72-3000
Small Epistaxis Balloon Catheter	Box of 1	72-3001

NETWORK NASAL SPLINTS


- Internal and external nasal splints.
- Airway splints. Silicone provides septal support, airway management and easier breathing.



EXTERNAL NASAL SPLINT

Supplied Non-Sterile

Description	Size	Pack	Code
External Nasal Splint, Malleable	Small	10	72-3020
External Nasal Splint, Malleable	Medium	10	72-3021
External Nasal Splint, Malleable	Large	10	72-3022



EXTERNAL NASAL SPLINT

Supplied Non-Sterile

Description	Size	Pack	Code
External Nasal Splint, Thermoplastic	Small	10	72-3025
External Nasal Splint, Thermoplastic	Medium	10	72-3026
External Nasal Splint - Custom Cut	Large	10	72-3027



AIRWAY SPLINT

Supplied Sterile

Description	Size	Pack	Code
Airway Splint	Standard	10 (Pairs)	72-3030



BIVALVE SPLINT

Supplied Sterile

Description	Size	Pack	Code
Bivalve Splint, Thin	Standard	10	72-3032
Bivalve Splint, Thick	Standard	10	72-3033
Bivalve Splint, Thin	Large	10	72-3034
Bivalve Splint, Thick	Large	10	72-3035




Description	Size	Pack	Code
Septal Button	32mm	Single	72-3040

Supplied Non-Sterile



VENTILATION TUBES

Dimensions key: **A** = Internal Diameter **B** = Inner Flange Diameter



Supplied Sterile

	Description	Material	A x B x C	Pack	Code
	T-Tube (12mm)	Silicone	1.1 x 9.0 x 12.0 mm	10 per box	74-100
	T-Tube (9mm)	Silicone	1.1 x 9.0 x 9.0 mm	10 per box	74-101
	T-Tube (6mm)	Silicone	1.1 x 9.0 x 6.0 mm	10 per box	74-102
	Paparella Tube II	Silicone	1.5 x 4.5 x 1.1 mm	10 per box	74-103
	Paparella Tube I	Silicone	1.0 x 2.1 x 1.1 mm	10 per box	74-104
	Shepard Grommets with wire	Fluoroplastic	1.1 x 2.3 x 1.5 mm	10 per box	74-106
	Shepard Grommets with wire	Fluoroplastic	1.1 x 2.3 x 1.5 mm	2 per pouch; 20 pouches	74-107
	Donaldson	Silicone	1.1 x 2.4 x 0.8 mm	10 per box	74-108
	Donaldson with tab	Silicone	1.1 x 2.4 x 0.8 mm	10 per box	74-109
	Shepard Grommets without wire	Fluoroplastic	1.1 x 2.3 x 1.5 mm	10 per box	74-110
	Shah without wire	Fluoroplastic	1.1 x 3.1 x 1.5 mm	10 per box	74-111
	Shah with wire	Fluoroplastic	1.1 x 3.1 x 1.5 mm	10 per box	74-112
	Straight Tube	Fluoroplastic	0.9 x 2.2 x 7.0 mm	10 per box	74-113
	Straight Tube	Fluoroplastic	1.1 x 2.2 x 7.0 mm	10 per box	74-114
	Armstrong	Fluoroplastic	1.1 x 2.6 x 4.0 mm	10 per box	74-115
	Armstrong (plain end)	Silicone	1.1 x 2.6 x 7.0 mm	10 per box	74-116
	Shepard with tail	Fluoroplastic	1.1 x 2.4 x 1.5 mm	10 per box	74-117
	Donaldson	Fluoroplastic	1.1 x 2.3 x 0.8 mm	10 per box	74-120
	Reuter Bobbin	Fluoroplastic	1.0 x 2.5 x 0.9 mm	10 per box	74-121
	Collar Button	Fluoroplastic	1.2 x 2.9 x 1.3 mm	10 per box	74-122
	Collar Button	Fluoroplastic	1.0 x 2.5 x 1.2 mm	10 per box	74-123
	Collar Button	Titanium	1.2 x 2.5 x 1.1 mm	10 per box	74-125

VARI-TUBES






70mm ANGLED VARI-TUBES

Supplied Sterile

	Dia. / Gauge	Hub Colour	Pack	Code
	2.0mm / 14G	Green	50 per box	74-1230
	1.6mm / 16G	White	50 per box	74-1231
	1.2mm / 18G	Pink	50 per box	74-1232
	0.9mm / 20G	Yellow	50 per box	74-1233
	0.7mm / 22G	Black	50 per box	74-1234
	0.5mm / 24G	Purple	50 per box	74-1236


70mm STRAIGHT VARI-TUBES

Supplied Sterile

	Dia. / Gauge	Hub Colour	Pack	Code
	2.0mm / 14G	Green	50 per box	74-1240
	1.6mm / 16G	White	50 per box	74-1241
	1.2mm / 18G	Pink	50 per box	74-1242
	0.9mm / 20G	Yellow	50 per box	74-1243
	0.7mm / 22G	Black	50 per box	74-1244
	0.5mm / 24G	Purple	50 per box	74-1246

80mm ANGLED VARI-TUBES

Supplied Sterile

	Dia. / Gauge	Hub Colour	Pack	Code
	2.0mm / 14G	Green	50 per box	74-1260

HOUSE ADAPTER WITH SUCTION CONTROL AND LUER CONNECTION

Supplied Sterile

	Pack	Code
	50 per box	74-2001

VARI-TUBES

14 GAUGE OINTMENT APPLICATOR

Supplied Sterile



Dia. / Gauge

Pack

Code

2.0mm / 14G

50 per box

74-1500

VARIFLEX SUCTION SET

Supplied Sterile



Description

Dia. / Gauge

Pack

Code

Variflex suction set complete with Zoellner type 2mm x 85mm, 10 degree shallow distal bend without suction control hole

2.0mm / 6Fr

50 per box

74-2004



Description

Dia. / Gauge

Pack

Code

Variflex suction set complete with Zoellner type 2mm x 85mm, 30 degree acute distal bend without suction control hole

2.0mm / 6Fr

50 per box

74-2005



Description

Dia. / Gauge

Pack

Code

Variflex suction set complete with Belucci type 2mm x 85mm proximal bend with suction control hole

2.0mm / 6Fr

50 per box

74-2006

FINE ENDS

Supplied Sterile









	Description	Hub Colour	Pack	Code
	Zoellner Fine End 18G, 1.0mm	Purple	50 per box	74-2020
	Zoellner Fine End 20G, 0.8mm	Orange	50 per box	74-2021
	Zoellner Fine End 22G, 0.65mm	Blue	50 per box	74-2022
	Zoellner Fine End 24G, 0.5mm	Green	50 per box	74-2023
	Zoellner Fine End 26G, 0.4mm	Yellow	50 per box	74-2024



VARI-TUBES

VARI-TUBE SUCTION HANDLES

Supplied Sterile

Description	Dia. / Gauge	Pack	Code
 <p>Brackmann type 1.6mm x 115mm, 40 degree bend, with suction control hole. The Brackmann type suction cannula has been developed to avoid trauma on delicate structures through strong suction. It is particularly useful for gentle suction around the facial nerve and labyrinthine fistulae and for removal of debris during Stapes, Myringoplasty and Ossiculoplasty surgery.</p>	1.6mm / 5Fr	10 per box	74-2008
 <p>Zoellner type 2mm x 85mm, 10° shallow distal bend without suction control hole</p>	2.0mm / 6Fr	50 per box	74-2011
 <p>Zoellner type 2mm x 85mm, distal bend acute 30° angle without suction control hole</p>	2.0mm / 6Fr	50 per box	74-2012
 <p>Zoellner type 2mm x 85mm, 10° shallow distal bend with suction control hole</p>	2.0mm / 6Fr	50 per box	74-2013
 <p>Zoellner type 2mm x 85mm, distal bend acute 30° angle with suction control hole</p>	2.0mm / 6Fr	50 per box	74-2014
 <p>Belucci type 2mm x 85mm, proximal bend with suction control hole</p>	2.0mm / 6Fr	50 per box	74-2015
 <p>Frazier Lampert Suction Tube 2mm x 125mm, proximal bend with suction control hole</p>	2.0mm / 6Fr	50 per box	74-2017
 <p>Frazier Lampert Suction Tube 3mm x 125mm, proximal bend with suction control hole</p>	3.0mm / 9Fr	50 per box	74-2050



NETWORK SKINMARKERS

- Skin-marking pens
- Supplied sterile

Supplied Sterile

Pack	Code
Box of 10	79-1005



NETBLADES

- Miniature blades
- Myringotomy blades
- Handles

MINIATURE BLADES

Supplied Sterile

	Description	Pack	Code
	Chisel Blade	12 per pack	80-162
	Round Tip	12 per pack	80-164
	Lance Tip	12 per pack	80-165
	60° Angled Tip	6 per pack	80-166
	Curved Tip	12 per pack	80-167
	Bent 60°, sharp all round Suitable for canal procedures	6 per pack	80-172

MYRINGOTOMY BLADES

Supplied Sterile


	Description	Pack	Code
	Myringotomy Blade, Downcutting	6 per pack	80-271
	Myringotomy Blade, Upcutting	6 per pack	80-272
	Myringotomy Blade, Spear Tip	6 per pack	80-273
	Myringotomy Blade, Lance Tip	6 per pack	80-274

HANDLES

Supplied Non-Sterile

	Description	Size	Pack	Code
	Round, Knurled handle	10 cm	single item	80-203K
	Round, Knurled handle	15.5 cm	single item	80-203KL
	Hexagonal handle	15.5 cm	single item	80-203HL

NOTES



Network Medical Products, The Coach House,
Phoenix Business Centre, Low Mill Road, Ripon,
North Yorkshire HG4 1NS, England

Tel: +44 (0) 1765 609555

Fax: +44 (0) 1765 608476

Email: marketing@networkmedical.co.uk

www.networkmedical.co.uk

