




DCR and Lachrymal Intubation

Supplied sterile for single use only

Sterilisation by gamma irradiation

Declared 5-year shelf life



Crawford-type Intubation Set



Probe length
150mm

Silicone tube length
320mm

1 per box

Code 51-918

Retrieval device for Crawford-type Intubation Set




Overall length
125mm

1 per box

Code 51-919

Canaliculus Intubation Set



Probe length
185mm

Silicone tube length
300mm

1 per box

Code 51-920

DCR Bodkin Straight



Bodkin length 45mm
20 gauge

Silicone rod length
400mm

Box of 3

Code 51-922

DCR Bodkin Angled



Bodkin length 45mm
20 gauge

Silicone rod length
400mm

Box of 3

Code 51-923

Network Medical Products Ltd.

The Coach House, Phoenix Business Centre,
Low Mill Road, Ripon, North Yorkshire HG4 1NS, UK.

Telephone: +44 1765 609555

Fax: +44 1765 608476

email: marketing@networkmedical.co.uk

www.networkmedical.co.uk

