




NETWORK

MEDICAL PRODUCTS



NETWORK ENT PRODUCTS

ENT Single Use Disposable Products

NETCELL

ENT PVA Sponge Products

NETCELL ENT PVA Sponge Products

| | |
|------------------------|-----|
| Nasal Packs | 3-5 |
| Sinus Packs..... | 6 |
| Ear Packs | 6 |
| Microspheres | 7 |
| Instrument Wipes | 7 |
| Headbands..... | 7 |

NETWORK ENT PRODUCTS ENT Single Use Disposable Products

| | |
|---|-------|
| Stapes Piston..... | 8 |
| Epistaxis Balloon Catheter | 8 |
| Nasal Splints | 9 |
| Ventilation Tubes | 10 |
| Aspiration Cannulas | 11-12 |
| Suction Tube Sets | 13 |
| Suction Tubes | 14-15 |
| Fine Ends..... | 15 |
| Oesophageal/Laryngeal Suction Tubes | 16-18 |
| Specialist Suction Tubes..... | 19 |
| Yankauer Suction Tubes | 20 |
| Single Use Disposable Instruments..... | 21-22 |

Network Medical has over 20 years experience as a fully integrated, designer and manufacture of high quality ENT surgical products.

Our reputation has been built on consistent and reliable product quality, that surgeons trust time and time again.

All of our products are manufactured in the UK, and have been designed to work in harmony with surgeons preferred techniques, without compromising on patient comfort.

“At a time when many companies in the medical sector have outsourced their production, Network Medical has stayed true to its core values, by retaining full start to end manufacturing control, at our UK based factory.”

Phillip Weston, CEO





NETCELL NASAL PACK

- For use after septal, turbinate and rhinoplastic surgery.
- Provides gentle, equal pressure within the nasal cavity when absorbing post-operative bleeding.
- Available with an airway tube which will enhance patient comfort.
- Slimpack developed for use in conjunction with a nasal splint following septal surgery.
- Layered slimpack permits the user to customise the thickness of the pack.



NETCELL NASAL PACK WITH TUBE

The Netcell Nasal Packing can be supplied with an integral airway tube. The airway tube benefits the patient through easier breathing and improved patient comfort, it improves venting and reduces the risk of obstruction.

Supplied Sterile

| Product | Dimensions L x Th x H (cm) | Pack | Code |
|---------------------------------|----------------------------|-----------|---------|
| Nasal Pack 3.0 cm | 3.0 x 1.5 x 2.0 cm | Box of 20 | 10-095 |
| Nasal Pack 4.5 cm | 4.5 x 1.5 x 2.0 cm | Box of 20 | 10-100 |
| Nasal Pack with Tube 4.5 cm | 4.5 x 1.5 x 2.0 cm | Box of 20 | 10-100T |
| Nasal Slimpack 4.5 cm | 4.5 x 1.0 x 2.0 cm | Box of 20 | 10-101 |
| Nasal Pack 6.0 cm | 6.0 x 1.5 x 2.0 cm | Box of 20 | 10-105 |
| Nasal Pack with Tube 6.0 cm | 6.0 x 1.5 x 2.0 cm | Box of 20 | 10-105T |
| Nasal Slimpack 6.0 cm | 6.0 x 1.0 x 2.0 cm | Box of 20 | 10-106 |
| Nasal Pack 8.0 cm | 8.0 x 1.5 x 2.0 cm | Box of 20 | 10-110 |
| Nasal Pack with Tube 8.0 cm | 8.0 x 1.5 x 2.0 cm | Box of 20 | 10-110T |
| Nasal Slimpack 8.0 cm | 8.0 x 1.0 x 2.0 cm | Box of 20 | 10-111 |
| Nasal Slimpack layered 8.0cm | 8.0 x 1.0 x 2.0 cm | Box of 20 | 10-111L |
| Nasal Slimpack with Tube 8.0 cm | 8.0 x 1.0 x 2.0 cm | Box of 20 | 10-111T |
| Nasal Pack 10.0 cm | 10.0 x 1.5 x 2.0 cm | Box of 20 | 10-115 |
| Nasal Slimpack 10.0cm | 10.0 x 1.0 x 2.0 cm | Box of 20 | 10-116 |




NETCELL ANATOMICAL PACK

- Clinically designed to fill the nasal cavity - providing gentle and equal support and allowing the efficient absorption of post operative bleeding.
- Available with an airway tube which will enhance patient comfort.

Supplied Sterile

| Product | Dimensions L x Th x H | Pack | Code |
|--------------------------------------|-----------------------|-----------|---------|
| Anatomical Pack 8.0 cm | 8.0 x 1.5 x 3.0 cm | Box of 10 | 10-120 |
| Anatomical Pack with Tube 8.0 cm | 8.0 x 1.5 x 3.0 cm | Box of 10 | 10-120T |
| Anatomical Slimpack 8.0 cm | 8.0 x 1.0 x 3.0 cm | Box of 10 | 10-125 |
| Anatomical Slimpack with Tube 8.0 cm | 8.0 x 1.0 x 3.0 cm | Box of 10 | 10-125T |



NETCELL EPISTAXIS PACK

- Immediate, effective control of bleeding.
- Fluid absorption causes the packing to expand producing a tamponading effect.
- Long Epistaxis for the simple quick packing of the entire nasal cavity, the most efficient way to contain and restrict blood loss.
- Short Epistaxis pack for the quick convenient treatment of anterior Epistaxis.
- Layered pack allows the packing to be adjusted to the individual patient.

Supplied Sterile

| Product | Dimensions L x Th x H | Pack | Code |
|---|-----------------------|-----------|---------|
| Epistaxis Pack Long 10.0 cm | 10.0 x 1.5 x 2.5 cm | Box of 10 | 10-140 |
| Epistaxis Slimpack Long 10.0 cm | 10.0 x 1.0 x 2.5 cm | Box of 10 | 10-141 |
| Epistaxis Slimpack Long Layered 10.0 cm | 10.0 x 1.0 x 2.5 cm | Box of 10 | 10-141L |
| Epistaxis Pack Medium 8.0 cm | 8.0 x 1.5 x 2.5 cm | Box of 10 | 10-145 |
| Epistaxis Slimpack Medium 8.0 cm | 8.0 x 1.0 x 2.5 cm | Box of 10 | 10-146 |
| Epistaxis Pack Short 5.5 cm | 5.5 x 1.5 x 2.5 cm | Box of 20 | 10-150 |

**SERIES 5000
SINGLE HYDRATION NASAL PACK**



SERIES 5000 combines a unique non-stick outer surface with Netcell packing which is manufactured from superior quality, ultra-fine pore, highly absorbant PVA material. This combination produces the most effective user friendly nasal packing with the benefits of single hydration.

- The soft covered surface provides a non-stick interface against tissue ingrowth.
- Once inserted and hydrated the pack will remain moist for up to 7 days without re-hydration.
- Pain and discomfort during removal is radically improved.
- Double compressed for easy insertion – more than 50% narrower than competitive products.
- The outer covering will accept water based or oil based medication.
- Will not pre-expand during insertion.
- Strung for convenience and supplied sterile.

“We believe this is the best PVA-ENT product in the market. Recent Studies have confirmed that by avoiding all contact with tissue results in the least trauma during removal. Patient comfort is a key attribute for this product.”

Supplied Sterile

| Product | Dimensions L x Th x H | Pack | Code |
|--|-----------------------|-----------|--------|
| Single Hydration Nasal Pack 6.0 cm | 6.0 x 1.5 x 2.0 cm | Box of 10 | 15-105 |
| Single Hydration Nasal Pack 8.0 cm | 8.0 x 1.5 x 2.0 cm | Box of 10 | 15-110 |
| Single Hydration Nasal Slimpack 8.0 cm | 8.0 x 1.0 x 2.0 cm | Box of 10 | 15-111 |



NETCELL SINUS PACK

- Developed as a post-operative packing for F.E.S.S.
- Designed to prevent lateralization of the middle turbinates post-operatively.
- Can be used to pack the Osteomeatal complex to leave the nasal passage clear for easier breathing.

Supplied Sterile

| Product | Dimensions L x Th x H | Pack | Code |
|----------------------------|-----------------------|---------------------------|--------|
| Sinus Pack Standard 2.5 cm | 2.5 x 1.20 x 2.0 cm | Pack of 2/box of 10 Packs | 20-200 |
| Sinus Pack Standard 3.5 cm | 3.5 x 1.20 x 1.2 cm | Pack of 2/box of 10 Packs | 20-210 |
| Sinus Pack Medium 3.5 cm | 3.5 x 0.90 x 1.2 cm | Pack of 2/box of 10 Packs | 20-211 |
| Sinus Pack Small 3.5 cm | 3.5 x 0.60 x 1.2 cm | Pack of 2/box of 10 Packs | 20-212 |



NETCELL EAR-WICK

- Ideal for delivering antibiotics when treating Otitis Externa.
- Instant patient relief when hydrated with medicated otic solution.
- Available as fenestrated to enable aeration

Supplied Sterile

| Product | Dimensions | Pack | Code |
|----------------------------|------------------|------------------------------|--------|
| Ear Wick 15 mm | 9 mm dia x 15 mm | 1 per pack / box of 50 packs | 30-300 |
| Ear Wick 12 mm | 7 mm dia x 12 mm | 1 per pack / box of 10 packs | 30-301 |
| Ear Wick Fenestrated 15 mm | 9 mm dia x 15 mm | 1 per pack / box of 50 packs | 30-305 |
| Ear Wick 24 mm | 9 mm dia x 24 mm | 1 per pack / box of 10 packs | 30-310 |
| Ear Wick Fenestrated 24 mm | 9 mm dia x 24 mm | 1 per pack / box of 10 packs | 30-315 |



NETCELL POST-OP EAR PACK

- Designed to prevent stenosis of the canal following surgery.
- Will not expand lengthwise and is suitable to support a graft following tympanoplasty.

Supplied Sterile

| Product | Dimensions | Pack | Code |
|-------------------------------|-------------------|------------------------------|--------|
| Ear Pack (Strung) small 15 mm | 9 mm dia x 15 mm | 1 per pack / box of 10 packs | 30-320 |
| Ear Pack without string 24 mm | 12 mm dia x 24 mm | 1 per pack / box of 10 packs | 30-340 |




NETCELL MICROSpheres

- Developed to control haemostasis in any micro surgical technique.
- The microspheres are hydrated with the appropriate medication and applied to the areas of bleeding.
- Middle ear procedures to remove particles and control haemostasis.

Supplied Sterile

| Product | Dimensions | Pack | Code |
|-------------------|------------|------------------------------|--------|
| Microsphere Small | 3 mm dia | Pack of 10 / box of 10 packs | 30-380 |
| Microsphere Large | 5 mm dia | Pack of 10 / box of 10 packs | 30-382 |



INSTRUMENT WIPES

Micro-surgical instruments are expensive and easy to damage. Our Instrument Wipes have been especially developed for the maintenance and care of these delicate surgical products. The wipes are soft but very strong, they will not snag and will easily remove the protein that weaken and affect the surface and edges of instruments. The life of a surgical instrument will be extended if it is cleaned thoroughly immediately after use. The wipes are supplied sterile for use during surgery.

| Product | Sterile / Non Sterile | Pack | Code |
|-----------------|-----------------------|---------------------|--------|
| Instrument Wipe | Sterile | Pack of 1/box of 20 | 35-900 |




HEADBANDS

Beneficial for surgeons operating in demanding situations, when moistened with water the headband will have an immediate cooling effect through evaporation.

Reduces the risk of contamination to the operating site by controlling and limiting perspiration.

| Product | Sterile / Non Sterile | Pack | Code |
|----------|-----------------------|---------------------|--------|
| Headband | Non Sterile | Pack of 1/box of 10 | 35-950 |



STAPES PISTON

Trim-to length Teflon (fluoroplastic) stapes piston

- Fluoroplastic shaft may be trimmed to length easily using a sharp knife.
- Loop can be opened and the plastic memory of the fluoroplastic provides secure closure.
- Supplied sterile in single packs.

Supplied Sterile

| Shaft Diameter | Loop Diameter | Length | Pack | Code |
|----------------|---------------|--------|-------------|--------|
| 0.4 mm | 0.6 mm | 6.0 mm | Single Pack | 71-400 |
| 0.4 mm | 0.6 mm | 5.0 mm | Single Pack | 71-401 |
| 0.4 mm | 0.6 mm | 4.5 mm | Single Pack | 71-402 |
| 0.4 mm | 0.6 mm | 4.0 mm | Single Pack | 71-403 |
| 0.4 mm | 0.6 mm | 3.5 mm | Single Pack | 71-404 |
| 0.6 mm | 0.6 mm | 6.0 mm | Single Pack | 71-405 |
| 0.6 mm | 0.6 mm | 5.0 mm | Single Pack | 71-406 |
| 0.6 mm | 0.6 mm | 4.5 mm | Single Pack | 71-407 |
| 0.6 mm | 0.6 mm | 4.0 mm | Single Pack | 71-408 |
| 0.6 mm | 0.6 mm | 3.5 mm | Single Pack | 71-409 |
| 0.8 mm | 0.6 mm | 6.0 mm | Single Pack | 71-410 |
| 0.8 mm | 0.6 mm | 5.0 mm | Single Pack | 71-411 |
| 0.8 mm | 0.6 mm | 4.5 mm | Single Pack | 71-412 |
| 0.8 mm | 0.6 mm | 4.0 mm | Single Pack | 71-413 |
| 0.8 mm | 0.6 mm | 3.5 mm | Single Pack | 71-414 |



EPISTAXIS BALLOON CATHETER

- The most effective control of posterior epistaxis.
- Enhanced control of Intra Nasal Bleeding with the dual balloon system.
- Easy placement assisted by the bevelled posterior rim design.
- Airway management and patient breathing much assisted through the integral airway tube.
- Patient comfort is helped by using soft, medical grade silicone.
- Inflatable balloon cuffs conform to nasal cavity.
- Balloon size: 30cc and 10cc

Supplied Sterile

| Description | Length | Pack | Code |
|----------------------------------|--------|----------|---------|
| Epistaxis Balloon Catheter Large | 120 mm | Box of 1 | 72-3000 |
| Epistaxis Balloon Catheter Small | 97 mm | Box of 1 | 72-3001 |

NETWORK NASAL SPLINTS


- Internal and external nasal splints.
- Airway splints. Provides septal support, airway management and easier breathing.



EXTERNAL NASAL SPLINT: MALLEABLE

Supplied Non-Sterile

| Description | Pack | Code |
|--|------|---------|
| External Nasal Splint, Malleable, Small | 10 | 72-3020 |
| External Nasal Splint, Malleable, Medium | 10 | 72-3021 |
| External Nasal Splint, Malleable, Large | 10 | 72-3022 |



EXTERNAL NASAL SPLINT: THERMOPLASTIC

Supplied Non-Sterile

| Description | Pack | Code |
|--|------|---------|
| External Nasal Splint, Thermoplastic, Small | 10 | 72-3025 |
| External Nasal Splint, Thermoplastic, Medium | 10 | 72-3026 |
| External Nasal Splint, Thermoplastic – Custom Cut, Large | 10 | 72-3027 |



AIRWAY SPLINT: SILICONE

Supplied Sterile

| Description | Pack | Code |
|-------------------------------------|------------|---------|
| Airway Splint – Pair left and right | 10 (Pairs) | 72-3030 |



BIVALVE SPLINT: FLUOROPLASTIC

Supplied Sterile

| Description | Pack | Code |
|--------------------------------|------|---------|
| Bivalve Splint, Standard Thin | 10 | 72-3032 |
| Bivalve Splint, Standard Thick | 10 | 72-3033 |
| Bivalve Splint, Large Thin | 10 | 72-3034 |
| Bivalve Splint, Large Thick | 10 | 72-3035 |



SEPTAL BUTTON: SILICONE


Supplied Sterile

| Description | Pack | Code |
|------------------------------|--------|---------|
| Septal Button, Diameter 32mm | Single | 72-3040 |


















VENTILATION TUBES

Dimensions key: **A** = Internal Diameter **B** = Inner Flange Diameter





Supplied Sterile

| Description | Material | A x B x C | Pack | Code |
|---|---------------|---------------------|----------------------------|--------|
|  T-Tube (12mm) | Silicone | 1.1 x 9.0 x 12.0 mm | 10 per box | 74-100 |
|  T-Tube (9mm) | Silicone | 1.1 x 9.0 x 9.0 mm | 10 per box | 74-101 |
|  T-Tube (6mm) | Silicone | 1.1 x 9.0 x 6.0 mm | 10 per box | 74-102 |
|  Paparella Tube II | Silicone | 1.5 x 4.5 x 1.1 mm | 10 per box | 74-103 |
|  Paparella Tube I | Silicone | 1.0 x 2.1 x 1.1 mm | 10 per box | 74-104 |
|  Shepard Grommets with wire | Fluoroplastic | 1.1 x 2.3 x 1.5 mm | 10 per box | 74-106 |
|  Shepard Grommets with wire | Fluoroplastic | 1.1 x 2.3 x 1.5 mm | 2 per pouch; 20 pouches | 74-107 |
|  Donaldson | Silicone | 1.1 x 2.4 x 0.8 mm | 10 per box | 74-108 |
|  Donaldson with tab | Silicone | 1.1 x 2.4 x 0.8 mm | 10 per box | 74-109 |
|  Shepard Grommets without wire | Fluoroplastic | 1.1 x 2.3 x 1.5 mm | 10 per box | 74-110 |
|  Shah without wire | Fluoroplastic | 1.1 x 3.1 x 1.5 mm | 10 per box | 74-111 |
|  Shah with wire | Fluoroplastic | 1.1 x 3.1 x 1.5 mm | 10 per box | 74-112 |
|  Armstrong | Fluoroplastic | 1.1 x 2.6 x 4.0 mm | 10 per box | 74-115 |
|  Shepard with tail | Fluoroplastic | 1.1 x 2.4 x 1.5 mm | 10 per box | 74-117 |
|  Donaldson | Fluoroplastic | 1.1 x 2.3 x 0.8 mm | 10 per box | 74-120 |
|  Reuter Bobbin | Fluoroplastic | 1.0 x 2.5 x 0.9 mm | 10 per box | 74-121 |
|  Collar Button | Fluoroplastic | 1.2 x 2.9 x 1.3 mm | 10 per box | 74-122 |
|  Collar Button | Fluoroplastic | 1.0 x 2.5 x 1.2 mm | 10 per box | 74-123 |

ASPIRATION TUBES

70mm ANGLED ASPIRATION CANNULA

Supplied Sterile

| | Dia. / Gauge | Hub Colour | Pack | Code |
|---|--------------|------------|------------|----------|
|  | 2.0mm / 14G | Green | 50 per box | 74-1230 |
|  | 1.7mm / 15G | Grey | 50 per box | 74-1230A |
|  | 1.6mm / 16G | White | 50 per box | 74-1231 |
|  | 1.4mm / 17G | Red | 50 per box | 74-1231A |
|  | 1.2mm / 18G | Pink | 50 per box | 74-1232 |
|  | 1.0mm / 19G | Cream | 50 per box | 74-1232A |
|  | 0.9mm / 20G | Yellow | 50 per box | 74-1233 |
|  | 0.7mm / 22G | Black | 50 per box | 74-1234 |
|  | 0.55mm / 24G | Purple | 50 per box | 74-1236 |
|  | 0.4mm / 27G | Grey | 50 per box | 74-1237 |

70mm STRAIGHT ASPIRATION CANNULA

Supplied Sterile

| | Dia. / Gauge | Hub Colour | Pack | Code |
|---|--------------|------------|------------|---------|
|  | 2.0mm / 14G | Green | 50 per box | 74-1240 |
|  | 1.6mm / 16G | White | 50 per box | 74-1241 |
|  | 1.2mm / 18G | Pink | 50 per box | 74-1242 |
|  | 0.9mm / 20G | Yellow | 50 per box | 74-1243 |
|  | 0.7mm / 22G | Black | 50 per box | 74-1244 |
|  | 0.55mm / 24G | Purple | 50 per box | 74-1246 |

80mm ANGLED ASPIRATION CANNULA

Supplied Sterile

| | Dia. / Gauge | Hub Colour | Pack | Code |
|---|--------------|------------|------------|----------|
|  | 2.0mm / 14G | Green | 50 per box | 74-1260 |
|  | 0.7mm / 22G | Black | 50 per box | 74-1261A |
|  | 1.0mm / 19G | Cream | 50 per box | 74-1262A |
|  | 1.4mm / 17G | Red | 50 per box | 74-1264 |
|  | 1.7mm / 15G | Grey | 50 per box | 74-1265 |
|  | 2.5mm / 13G | Clear | 50 per box | 74-1267 |
|  | 3.0mm / 11G | Blue | 50 per box | 74-1268 |
|  | 4.0mm / 8G | Clear | 50 per box | 74-1270 |


100mm ANGLED ASPIRATION CANNULA


Supplied Sterile





| | Dia. / Gauge | Hub Colour | Pack | Code |
|---|--------------|------------|------------|---------|
|  | 1.4mm / 17G | Red | 50 per box | 74-1282 |
|  | 2.0mm / 14G | Green | 50 per box | 74-1290 |
|  | 2.5mm / 13G | Clear | 50 per box | 74-1294 |
|  | 3.0mm / 11G | Blue | 50 per box | 74-1297 |

HOUSE ADAPTER WITH SUCTION CONTROL

Supplied Sterile






| | Pack | Code |
|---|------------|---------|
|  | 50 per box | 74-2001 |

| | | | | Supplied Sterile |
|--|--------------|------------|---------|------------------|
| Description | Dia. / Gauge | Pack | Code | |
|  <p>Ointment Applicator supplied connected with a 2.5cc syringe</p> | 2.0mm / 14G | 50 per box | 74-1500 | |






| | | | | Supplied Sterile |
|---|--------------|------------|---------|------------------|
| SUCTION SET | | | | |
| Description | Dia. / Gauge | Pack | Code | |
|  <p>Suction set complete with house adapter with suction control hole</p> | | 50 per box | 74-2003 | |
| Description | Dia. / Gauge | Pack | Code | |
|  <p>Suction set with Zoellner type 2mm x 85mm, 10 degree shallow distal bend without suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2004 | |
| Description | Dia. / Gauge | Pack | Code | |
|  <p>Suction set complete with Zoellner type 2mm x 85mm, 30 degree acute distal bend without suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2005 | |
| Description | Dia. / Gauge | Pack | Code | |
|  <p>Suction set with Belucci type 2mm x 85mm proximal bend with suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2006 | |

SUCTION TUBES

Supplied Sterile






| Description | Dia. / Gauge | Pack | Code |
|---|--------------|------------|----------|
|  <p>Zoellner type 2mm x 85mm, 10° shallow distal bend without suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2011 |
|  <p>Zoellner type 2mm x 85mm, distal bend acute 30° angle without suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2012 |
|  <p>Zoellner Pen Grip type 2mm x 85mm, distal bend acute 30 deg angle without suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2012P |
|  <p>Zoellner type 2mm x 85mm, 10° shallow distal bend with suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2013 |
|  <p>Zoellner type 2mm x 85mm, distal bend acute 30° angle with suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2014 |

Supplied Sterile


| Description | Dia. / Gaug | Pack | Code |
|--|-------------|------------|-----------------|
|  <p>Zoellner Pen Grip type 2mm x 85mm, distal bend acute 30 deg angle with suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2014P |
|  <p>Belucci type 2mm x 85mm, proximal bend with suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2015 |
|  <p>Zoellner type 2mm x 100mm, medial bend acute 30 deg angle without suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2016 |
|  <p>Zoellner Pen Grip type 2mm x 100mm, medial bend acute 30 deg angle without suction control hole</p> | 2.0mm / 6Fr | 50 per box | 74-2016P |
|  <p>Magill Type 3mm x 110mm, 10 degree angle, shallow distal bend, with suction control hole</p> | 3.0mm / 9Fr | 50 per box | 74-2010 |

FINE ENDS

Supplied Sterile

| Description | Hub Colour | Pack | Code |
|---|------------|------------|----------------|
|  <p>Fine End 18G, 1.0mm</p> | Purple | 50 per box | 74-2020 |
|  <p>Fine End 20G, 0.8mm</p> | Orange | 50 per box | 74-2021 |
|  <p>Fine End 22G, 0.65mm</p> | Blue | 50 per box | 74-2022 |
|  <p>Fine End 24G, 0.55mm</p> | Green | 50 per box | 74-2023 |
|  <p>Fine End 26G, 0.4mm</p> | Yellow | 50 per box | 74-2024 |


OESOPHAGEAL/LARYNGEAL SUCTION TUBES



Supplied Sterile

| Description | Length | Dia. / Gauge | Pack | Code |
|---|--------|--------------|------------|----------------|
| Oesophageal / Laryngeal Suction Tubes. 2mm x 125mm, proximal bend with suction control hole | 125mm | 2.0mm / 6Fr | 50 per box | 74-2017 |
| Oesophageal / Laryngeal Suction Tubes. 2mm x 165mm Vari-Tube. 30 Deg angle, proximal bend. With Suction Control | 165mm | 2.0mm / 6Fr | 25 per box | 74-2029 |
| Oesophageal / Laryngeal Suction Tubes. 2mm x 215mm Vari-Tube. 30 Deg angle, proximal bend. With Suction Control | 215mm | 2.0mm / 6Fr | 25 per box | 74-2030 |
| Oesophageal / Laryngeal Suction Tubes. 2mm x 275mm Vari-Tube. 30 Deg angle, proximal bend. With Suction Control | 275mm | 2.0mm / 6Fr | 25 per box | 74-2033 |
| Oesophageal / Laryngeal Suction Tubes. 2mm x 350mm Vari-Tube. 30 Deg angle, proximal bend. With Suction Control | 350mm | 2.0mm / 6Fr | 25 per box | 74-2037 |


OESOPHAGEAL/LARYNGEAL SUCTION TUBES



Supplied Sterile

| Description | Length | Dia. / Gauge | Pack | Code |
|---|--------|--------------|-----------|-----------------|
| Frazier Lampert Type. 3mm x 125mm. Proximal bend. With Suction Control | 125mm | 3.0mm / 9Fr | Box of 50 | 74-2050 |
| Oesophageal / Laryngeal Suction Tubes. 3mm x 150mm. 30 Deg angle, proximal bend. With Suction Control | 150mm | 3.0mm / 9Fr | Box of 25 | 74-2050A |
| Oesophageal / Laryngeal Suction Tubes. 3mm x 165mm. 30 Deg angle, proximal bend. With Suction Control | 165mm | 3.0mm / 9Fr | Box of 25 | 74-2051 |
| Oesophageal / Laryngeal Suction Tubes. 3mm x 250mm. 30 Deg angle, proximal bend. With Suction Control | 250mm | 3.0mm / 9Fr | Box of 25 | 74-2052 |
| Oesophageal / Laryngeal Suction Tubes. 3mm x 275mm. 30 Deg angle, proximal bend. With Suction Control | 275mm | 3.0mm / 9Fr | Box of 25 | 74-2053 |
| Oesophageal / Laryngeal Suction Tubes. 3mm x 350mm. 30 Deg angle, proximal bend. With Suction Control | 350mm | 3.0mm / 9Fr | Box of 25 | 74-2057 |
| Oesophageal / Laryngeal Suction Tubes. 3mm x 380mm. 30 Deg angle, proximal bend. With Suction Control | 380mm | 3.0mm / 9Fr | Box of 25 | 74-2058 |

OESOPHAGEAL/LARYNGEAL SUCTION TUBES











Supplied Sterile

| Description | Length | Dia. / Gauge | Pack | Code |
|---|--------|--------------|-----------|----------------|
| Oesophageal / Laryngeal Suction Tubes. 4mm x 120mm. 30 Deg angle, proximal bend. With Suction Control | 120mm | 4.0mm / 12Fr | Box of 25 | 74-2077 |
| Oesophageal / Laryngeal Suction Tubes. 4mm x 150mm. 30 Deg angle, proximal bend. With Suction Control | 150mm | 4.0mm / 12Fr | Box of 25 | 74-2079 |
| Oesophageal / Laryngeal Suction Tubes. 4mm x 165mm. 30 Deg angle, proximal bend. With Suction Control | 165mm | 4.0mm / 12Fr | Box of 25 | 74-2080 |
| Oesophageal / Laryngeal Suction Tubes. 4mm x 125mm. 30 Deg angle, proximal bend. With Suction Control | 125mm | 4.0mm / 12Fr | Box of 25 | 74-2081 |
| Oesophageal / Laryngeal Suction Tubes. 4mm x 215mm. 30 Deg angle, proximal bend. With Suction Control | 215mm | 4.0mm / 12Fr | Box of 25 | 74-2082 |
| Oesophageal / Laryngeal Suction Tubes. 4mm x 275mm. 30 Deg angle, proximal bend. With Suction Control | 275mm | 4.0mm / 12Fr | Box of 25 | 74-2083 |
| Oesophageal / Laryngeal Suction Tubes. 4mm x 350mm. 30 Deg angle, proximal bend. With Suction Control | 350mm | 4.0mm / 12Fr | Box of 25 | 74-2087 |

SPECIALIST SUCTION TUBES



Supplied Sterile

| Description | | Dia. / Gauge | Pack | Code |
|--|---|--------------|------------|---------|
|  <p>Brackmann type with end cap 1.6mm x 115mm, 40 degree bend, with suction control hole.</p> |  | 1.6mm / 5Fr | 10 per box | 74-2008 |
|  <p>Brackmann type without end cap 1.6mm x 115mm, 40 degree bend, with suction control hole.</p> <p>Avoid trauma on delicate structures for gentle suction around the facial nerve and labyrinthine fistulae and for removal of debris during Stapes, Myringoplasty and Ossiculoplasty surgery.</p> |  | 1.6mm / 5Fr | 10 per box | 74-2009 |
|  <p>F.E.S.S sinus suction tube 3mm x 100mm, distal bend, 90 degree angle with olive tip and suction control hole</p> |  | 3.0mm / 9Fr | 25 per box | 74-2060 |
|  <p>F.E.S.S sinus suction tube 3mm x 100mm, distal bend, 65 degree angle with olive tip and suction control hole</p> |  | 3.0mm / 9Fr | 25 per box | 74-2063 |

YANKAUER SUCTION TUBES



Features

- Ergonomically designed to feel natural in the hand
- Choice of type, size and angle
- Manufactured with clear tips
- Ribbed universal connector
- Supplied with open tips

Benefits


- Gives superior balance and ease of use
- Offers suction handle suitable for all procedures
- Permits visual monitoring of suction fluids
- Fits securely on all types of tubing
- High volume suction achieved

Supplied Sterile

| Description | Dia. / Gauge | Pack | Code |
|--|--------------|-----------|---------|
| 8mm Open Tip, Small Curve with Suction Control Hole, with vents | 8mm | Box of 50 | 74-1211 |
| 8mm Open Tip, Small Curve without Suction Control Hole | 8mm | Box of 50 | 74-1212 |
| 6mm Flange Tip, Medium Curve with Suction Control Hole, with vents | 6mm | Box of 50 | 74-1213 |
| 6mm Flange Tip, Medium Curve without Suction Control Hole | 6mm | Box of 50 | 74-1214 |
| 10mm Bulbous Tip, Large Curve with Suction Control Hole, with vents | 10mm | Box of 50 | 74-1215 |
| 10mm Bulbous Tip, Large Curve without Suction Control Hole, with vents | 10mm | Box of 50 | 74-1216 |

ROSEN AURAL SPECULA

(UK and European Market Only)



Matt black finish for improved visibility under the microscope. Supplied sterile for single use ensuring optimal hygiene.

Sterile

| Diameter | Pack | Code |
|----------|------------|---------|
| 3mm ○ | 50 per box | 81-1000 |
| 4mm ○ | 50 per box | 81-1001 |
| 5mm ○ | 50 per box | 81-1002 |
| 6mm ○ | 50 per box | 81-1003 |
| 7mm ○ | 50 per box | 81-1004 |

TUMARKIN AURAL SPECULA

(UK and European Market Only)



Matt black finish for improved visibility under the microscope. Supplied sterile for single use ensuring optimal hygiene.

Sterile

| Diameter | Pack | Code |
|----------|------------|---------|
| 3.5mm ○ | 50 per box | 81-1008 |
| 4.5mm ○ | 50 per box | 81-1009 |
| 5.5mm ○ | 50 per box | 81-1010 |
| 6.0mm ○ | 50 per box | 81-1011 |
| 6.5mm ○ | 50 per box | 81-1012 |
| 7.0mm ○ | 50 per box | 81-1013 |



JOBSON HORNE PROBE

(UK and European Market Only)



Supplied Sterile

| Description | Pack | Code |
|---|------------|---------|
| Rounded surface with no sharp edges. Lightweight for precise control. | 50 per box | 81-1020 |



THUDICHUMS

(UK and European Market Only)

Smooth satin finish for no reflection. Positive spring pressure. No sharp edges. Lightweight for added patient comfort.

Supplied Sterile

| Size | Pack | Code |
|--------|------------|---------|
| Small | 20 per box | 81-1030 |
| Medium | 20 per box | 81-1031 |
| Large | 20 per box | 81-1032 |

CROCODILE FORCEPS

(UK and European Market Only)



Precision machined and individually assembled, with a positive balance and clean-cut serrated jaws that give a reliable grip and no snagging.

Supplied Sterile

| Description | Jaws | Pack | Code |
|--|------------|------------|----------------|
| Micro crocodile forceps, black finish | Extra fine | 20 per box | 81-2000 |
| Standard crocodile forceps, satin finish | Standard | 20 per box | 81-2012 |



VANNAS SCISSORS

(UK and European Market Only)

Suitable for micro-suture and Septal Rhinoplasty procedures. Stainless steel with straight blades, supplied sterile



Supplied Sterile

| Description | Pack | Code |
|--------------------------|------------|----------------|
| Vannas scissors straight | 20 per box | 81-3000 |



Network Medical Products Ltd.
Coronet House, Kearsley Road,
Ripon, North Yorkshire, HG4 2SG, UK.

Tel: +44 (0) 1765 609555 Fax: +44 (0) 1765 608476

email: info@networkmedical.co.uk

www.networkmedical.co.uk

Distributed by: