


The NETCELL™ Series 5000™ PVA Nasal Pack has been clinically proven to reduce trauma and improve patient comfort¹

¹Source International Rhinology 47, 41-44. 2009 A randomised controlled trial comparing Rapid Rhino Mannheim and NETCELL Series 5000 packs following routine surgery



NETCELL Series 5000 (NET)

Figure 3. Box plots of change in VAS pain score for each pack type. The p value for the difference between these 2 data sets is 0.008.

- Unique patent protected nonstick and soft outer surface prevents tissue in-growth
- Specifically designed for use after septal, turbinate and rhinoplasty surgery
- Provides gentle, equal pressure within the nasal cavity when stenting post-operative bleeding
- Remains hydrated for up to 7 days without re hydration
- Double compressed for easy insertion (50% slimmer than competitor products)
- Outer coating can be used in conjunction with water or oil based medication
- Supplied sterile, single use only, declared 5 year shelf life

Product Description	Dimension LXTHxH	Pack Size	Code
Single hydration nasal packs 6.0cm	6.0 x 1.5 x 2.0cm	Box of 10	15-105
Single hydration nasal pack 8.0cm	8.0 x 1.5 x 2.0cm	Box of 10	15-110
Single hydration nasal pack slimpack 8.0cm	6.0 x 1.0 x 2.0cm	Box of 10	15-111