



### Product Overview

The NETCELL® PVA Instrument Wipes have been specifically developed to effectively clean and extend the life of delicate micro surgical instruments.

- Quickly removes protein, blood and particles without damage to the instrument
- Manufactured from ultra-smooth, biocompatible PVA sponge that will not snag
- Soft, strong, flexible, fibre and snag free material
- Supplied sterile, single use only, declared 5 year shelf life

### PVA Instrument Wipe Size Options

Product Description	Format	Pack Size	Product Code
Instrument Wipe	Standard	20	35-900

### Material Specification

---

#### Product Component

#### Specification

Instrument Wipe	100% Polyvinyl Alcohol Sponge (PVA)
Pouch Packaging	Tyvek/Film
Outer Box/Carton Packaging	500 Micron Printed White Boxboard

---

### Intended Use

PVA instrument wipes are single use, disposable cleaning wipe, for the removal of blood, particles and protein from surgical instruments.

### Sterilisation

Products are sterilised by Gamma irradiation from a Cobalt 60 source in accordance with a validated 25 kGy cycle. Sterilisation is carried out in accordance with the requirements of ISO 11137-1 (current) and ISO 11137-2 (current) and the 25 kGy dose is substantiated by VD<sub>25</sub> Method Max testing.

### Instructions for Use

Instructions for use are not provided as this product can be used safely without.

### Conformity to the European Directives

Instrument wipes are not considered to be medical devices as they are not intended for patient contact. Sold for single use only, the sterilisation process is covered by the ISO 13485:2012 certification carried by Network Medical Products Ltd.