

Aspiration Cannula Range Specification Sheet



Product Overview

The NETWORK ENT® range of Aspiration Cannulas have been specifically designed to offer precision suction for intricate ENT procedures and are suitable for use with most suction systems and house adapters.

- Ultra- thin walled tubing made from high-grade Japanese steel, allows a smaller OD for increased precision and control.
- Ergonomically designed to provide superior balance and improved ease of use
- Available in a wide range of sizes in both straight and angled format
- Compatible with the NETWORK ENT House Adapter with suction control (74-2001)
- Supplied sterile, single use only, declared 5 year shelf life
- Manufactured in the UK

Aspiration Cannula Range Specification Sheet

70mm Angled Size Options

Product Size	Diameter/Gauge	Pack Size	Hub Colour	Product Code
70mm Angled	2.0mm/14G	50 per pack	Green	74-1230
70mm Angled	1.7mm/15G	50 per pack	Grey	74-1230A
70mm Angled	1.6mm/16G	50 per pack	White	74-1231
70mm Angled	1.4mm/17G	50 per pack	Red	74-1231A
70mm Angled	1.2mm/18G	50 per pack	Pink	74-1232
70mm Angled	1.0mm/19G	50 per pack	Cream	74-1232A
70mm Angled	0.9mm/20g	50 per pack	Yellow	74-1233
70mm Angled	0.7mm/22G	50 per pack	Black	74-1234
70mm Angled	0.55mm/24G	50 per pack	Purple	74-1236

80mm Angled Size Options

Product Size	Diameter/Gauge	Pack Size	Hub Colour	Product Code
80mm Angled	3.0mm/11G	50 per pack	Blue	74-1268
80mm Angled	2.5mm/13G	50 per pack	Clear	74-1267
80mm Angled	2.0mm/14G	50 per pack	Green	74-1260
80mm Angled	1.7mm/15G	50 per pack	Grey	74-1265
80mm Angled	1.4mm/17G	50 per pack	Red	74-1264
80mm Angled	1.0mm/19G	50 per pack	Cream	74-1262A
80mm Angled	0.7mm/22G	50 per pack	Black	74-1261A

Aspiration Cannula Range Specification Sheet

100mm Angled Size Options

Product Size	Diameter/Gauge	Pack Size	Hub Colour	Product Code
100mm Angled	3.0mm/11G	50 per pack	Blue	74-1297
100mm Angled	2.5mm/13G	50 per pack	Clear	74-1294
100mm Angled	2.0mm/14G	50 per pack	Green	74-1290

Material Specification

Product Component	Specification
Aspiration Cannula Tube	Stainless Steel
Aspiration Cannula Hub	Polypropylene (PP)
Pouch Packaging	Tyvek/Film
Outer Box/Carton Packaging	500 Micron White Boxboard

Intended Use

The NETWORK ENT® Aspiration Cannulas are designed specifically for ENT procedures and are intended for the removal of excess fluids and debris from the surgical site.

The NETWORK ENT® Aspiration Cannula can be directly attached to a House Adapter or other adapter handle with a push fit connection. Each Aspiration Cannula is designated for 'SINGLE USE ONLY' and should not be re-sterilised by any method or re-used.

Instructions for Use

The instructions for use and suggested surgical technique are supplied in the form of a multi-lingual leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are explained in each language. A leaflet is supplied with each product.

Sterilisation