










Beaver™ Otolaryngology Blades

Beaver™ is part of the BVI Group of companies






NETWORK
ENT PRODUCTS










BEAVER™ MINI-BLADES

	Description	Pack	Code
	Beaver® 6400 Mini-Blade® Grindless, rounded tip, sharp on one side, green	20 per box	Beaver6400
	Beaver® 6600 Mini-Blade® Grindless, lamellar, angled 60 degrees, smooth undersurface, sharp all around, bevel up, black	3 per box	Beaver6600
	Beaver® 6610 Mini-Blade® Grindless, lamellar, straight, smooth undersurface, sharp all around, bevel up, black	20 per box	Beaver6610
	Beaver® 6700 Mini-Blade® Grindless, smaller alternative to #15 general blade, orange	20 per box	Beaver6700
	Beaver® 6900 Mini-Blade® Grindless, straight, sharp all around, 180° cutting surface, double bevel blade, blue	20 per box	Beaver6900
	Beaver® 6910 Mini-Blade® Grindless, angled 60 degrees, sharp all around, 180° cutting surface, double bevel blade, blue	6 per box	Beaver6910

Directions for Assembly

<p>1</p>  <p>Hold the Beaver® Mini-Blade® tray between thumb and pointer finger.</p>	<p>2</p>  <p>Thread the Beaver® handle in a clockwise motion, onto the threaded tip.</p>	<p>3</p>  <p>Continue to thread handle until the thread is fully engaged and the Mini-Blade® is released from the tray clips.</p>	<p>4</p>  <p>Lift to remove Beaver® Mini-Blade® from tray.</p>
--	--	--	--

ENT/HEAD & NECK PROCEDURE BLADES

	Description	Pack	Code
	Knife Straight, sharp all around, 180° cutting surface double bevel blade	3 per box	374769
	Needle Blade 3.0mm cutting edge	3 per box	375910
	Needle Blade 5.0 mm cutting edge	3 per box	375920
	Needle Blade 7.0 mm cutting edge	3 per box	375930
	Needle Blade 1.5mm cutting edge	3 per box	375940
	Beaver Mini-Blade Sharp on two sides, pointed tip	6 per box	376300
	Beaver Mini-Blade Rounded narrow tip, sharp on one side	24 per box	376401
	Tonsillectomy Blade Sickle sharp	24 per box	376800
	Beaver Mini-Blade Angled 60°, sharp all around, 180° cutting surface double bevel blade	6 per box	376910

BEAVER™ HANDLES

Beaver Handles will accommodate all Beaver blades, supplied non-sterile.

	Description	Pack	Code
	Round, Knurled, Stainless 10.0 cm L x .50 cm D	1 per box	374310
	Round, Knurled, Stainless 13.0 cm L x .50 cm D	1 per box	371350

MYRINGOTOMY BLADES AND KNIVES

	Description	Pack	Code
	Myringotomy Blade Lance tip, downward cut, offset 45°	6 per box	377100
	Myringotomy Blade Lance tip, upward cut, offset 45°	6 per box	377110
	Myringotomy Blade Spear tip, offset 45°, narrow shaft	6 per box	377120
	Myringotomy Blade Lance tip, downward cut, narrow shaft, offset 45°	6 per box	377121
	Myringotomy Knife Spear tip, offset 45°, narrow shaft	3 per box	377122
	Myringotomy Knife Lance tip, downward cut, narrow shaft, offset 45°	3 per box	377124

TYMPANOPLASTY/STAPEDECTOMY BLADES

	Description	Pack	Code
	Tympoplasty Blade 60° angled, sharp all around 2.5mm W, bevel down	3 per box	377200
	Tympoplasty Blade Straight, sharp all around 2.5mm W	3 per box	377210
	Sickle Edge Blade Straight 84mm L 5.0mm cutting surface	3 per box	377300
	Tympoplasty Blade 50° angled, sharp all around 2.0mm W, bevel down	3 per box	374451

All Beaver blades are supplied sterile



Network Medical Products Ltd.

Coronet House, Kearsley Road,
Ripon, North Yorkshire, HG4 2SG, UK.

Tel: +44 (0) 1765 609555 Fax: +44 (0) 1765 608476

email: info@networkmedical.co.uk

www.networkmedical.co.uk

Distributed by:

All the products in this catalogue are manufactured in accordance with the requirements of ISO13485:2003 and the Medical Devices Directive 93/42/EEC and, where appropriate, are authorised to bear the CE mark.