



Oesophageal/Laryngeal Suction Tubes

The NETWORK ENT™ Single Use Suction Tube range has been specifically designed to provide precision and control during ENT procedures.

- Available in 2mm, 3mm and 4mm up to 380mm in length
- Manufactured with consistent lumen through the bend to reduce the risk of blockages
- Constructed from high quality Japanese steel for superior thin walled tubing
- Ergonomically balanced and lightweight handle for unrivaled maneuverability and control
- Supplied sterile, single use only, declared 5 year shelf life
- Manufactured in the UK

