

CORONET

DMEK Guarded Trephine Punch

More than just an excellent cut...

Now available
in 10.00mm



PATENT PROTECTED

 **NETWORK**
MEDICAL PRODUCTS | AN
INNOVIA
MEDICAL
COMPANY

CORONET

The CORONET® DMEK Guarded Trephine Punch offers more than just an excellent cut.

Designed to meet specific demands of DMEK surgery the punch supports instrument position during tissue manipulation by providing an ergonomic and anatomically shaped platform.



CORONET® DMEK Guarded Trephine Punch 51-855

KEY FEATURES

- Sharp, clean cutting guarded trephine with exclusive CORONET® Cathedral blade technology
- Pre-set depth to score through Descemet's membrane with a simple press
- Patent protected visualisation windows
- Stable platform base
- Flexible non-restrictive vacuum tubing
- Biocompatible and non-glare blade
- Supplied sterile, single use, 5 year shelf life
- Supplied one per box
- Manufactured in the UK

TREPINE SIZE (mm)

Size	Code
7.50	51-855-7.50
8.00	51-855-8.00
8.50	51-855-8.50
9.00	51-855-9.00
9.50	51-855-9.50
10.00	51-855-10.00

Network Medical Products Ltd. | Coronet House | Kearsley Road | Ripon | North Yorkshire | HG4 2SG | UK
Tel. +44 1765 609555 | www.networkmedical.co.uk

USA Enquiries:
Innovia Medical | Summit Medical 815 Vikings Parkway | Suite 100, |St. Paul, MN | 55121 | USA
Tel+1 888 2292875 | 651 7893939
orders@innoviamedical.com

SI_COR_005 REV 04 04/19