


## Single Use Endonasal Sinus Suction Tubes

Specifically designed for entering through the natural sinus ostium during Functional Endoscopic Sinus Surgery (F.E.S.S)



- Olive shaped suction tip for atraumatic entrance into the sinus ostium
- Curved tip for maximum visibility and access to difficult areas
- Constructed from high quality Japanese steel for superior thin walled tubing
- Manufactured without bending or pinching to create a continuous radius, reducing blockages
- Ergonomically balanced and lightweight handle for improved maneuverability and control
- Supplied Sterile, Single Use Only, Declared 5 Year Shelf Life
- Manufactured in the UK

### Product Description

### Pack Size

### Code

F.E.S.S sinus suction tube, 3mm x 85mm, distal bend 80 degree angle with olive tip and suction control hole

25

74-2060

F.E.S.S sinus suction tube, 3mm x 85mm, distal bend 65 degree angle with olive tip and suction control hole

25

74-2063