

CORONET

DALK Clear-Vue™ Recipient Vacuum Trepine

Controlled DALK incision in three simple steps



PATENT PROTECTED

 **NETWORK**
MEDICAL PRODUCTS | AN INNOVIA
MEDICAL
COMPANY

CORONET

No need to calibrate or measure trephine depth*

Controlled DALK incision in three simple steps

1. Position trephine



2. Remove centration tube



3. Insert guarded trephine



*Fixed depth blade will cut once inserted

The CORONET® DALK Clear-Vue™ Recipient Vacuum Trephine 51-844

KEY FEATURES

- Centration tube with cross hair reduces parallax
- Guarded trephine blade with exclusive CORONET® technology
- Soft adaptive skirt conforms to individual eye topography
- Light weight body sits on limbus
- Secure vacuum formed away from incision
- 360° visualisation and ergonomic finger grips
- Supplied sterile, single use only, one per box, 5-year shelf life
- Manufactured in the UK

TREPHINE SIZE (mm)

Size	Code
7.50	51-844-7.50
8.00	51-844-8.00
8.50	51-844-8.50
9.00	51-844-9.00

Network Medical Products Ltd. | Coronet House | Kearsley Road | Ripon | North Yorkshire | HG4 2SG | UK
Tel. +44 1765 609555 | www.networkmedical.co.uk

USA Enquiries:
Innovia Medical | Summit Medical 815 Vikings Parkway | Suite 100, | St. Paul, MN | 55121 | USA
Tel+1 888 2292875 | 651 7893939
orders@innoviamedical.com

SI_COR_011 REV 01 04/19