




DCR Bodkins & Lacrimal Intubation

Manufactured in the UK





DCR Bodkins & Lacrimal Intubation

Crawford-type Intubation Set

Complete set for Dacryocystorhinostomy, supplied with retrieval device. Retrieval device also available separately.

KEY FEATURES

- Probe length 150mm
- Silicone tube length 320mm
- Retrieval device 105mm length



51-918

Canaliculus Intubation Set

KEY FEATURES

- Probe length 114mm
- Silicone tube length 300mm
- Supplied sterile, single use only, one per box.
- Declared 5 year shelf life



51-920

DCR Bodkins Straight and Angled

KEY FEATURES

- Bodkin length 45mm with 20 gauge diameter
- Silicone rod length 400mm
- Available in straight or angled format



51-922



51-923

All products are supplied sterile, single use with a 5 year declared shelf life.

Product	Pack	Code
Retrieval device for Crawford-Type Intubation	Box of 5	51-917
Crawford-Type Intubation Set	Box of 1	51-918
Canaliculus Intubation Set	Box of 1	51-920
DCR Bodkins Straight	Box of 3	51-922
DCR Bodkins Angled	Box of 3	51-923