

Reusable Innovia Face Visor Product Specification



The Design Checklist

- Manufactured to highest PPE regulation requirements**
 CE Category III 2016/425 as per Annex II (Notified Body Satra Technology Europe Ltd)
- Reusable design reduces environmental impact**
 Allows for up to 5 days of safe wear by the same individual
- N rated hydrophilic anti-fog coating**
 Prevents condensation forming for an optically clear screen and maximum visibility (compliant with EN166:2001 Clause 7.3.2)
- Permanent marker compatible Visor includes space for individual name and customisation
- Pre-printed 5-day tally tick box feature allows wearer to keep track of number of days used and when to dispose of and replace
- Ultra-soft grooved foam headrest for increased comfort and removal of dermal heat
- Designed to provide face protection from aerosols and minimise airborne cross-contamination
- Meets all applicable requirements of EN 166:2001 Personal Eye Protection
- Provides top, side, and front face protection for the healthcare professional
- 22.8cm long screen provides maximum protection
- Elastic headband provides secure retention and easy application
- Can be worn over prescription eyewear and facemasks
- Anti-glare performance for optimum clarity within clinical environment
- Lightweight and comfortable
- Latex-free
- Manufactured in an ISO Class 9 or higher UK cleanroom
- Non-sterile, individually wrapped



DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

DTR Medical Ltd, 17 Clarion Court, Enterprise Park, Swansea SA6 8RF, UK
 Tel +44(0) 1792 79 79 10 Fax +44(0) 1792 79 79 55
 marketing@dtrmedical.co.uk | www.dtrmedical.com



**NETWORK
MEDICAL PRODUCTS**

AN
INNOVIA
MEDICAL
COMPANY

Coronet House, Kearsley Road, Ripon, North Yorkshire, HG4 2SG, UK
 Tel: +44 (0) 1765 609555 Fax: +44 (0) 1765 608476
 info@networkmedical.co.uk | www.networkmedical.co.uk

Reusable Innovia Face Visor Product Specification

Intended Use

The device is intended to provide top, side, and front face protection from aerosols and to minimise airborne cross-contamination.

This protective Face Visor is manufactured for COVID-protection only. This protective Face Visor is intended to provide protection against liquid splashes.

This Visor is not resistant against spark, flame, or impact.

Instructions For Use (IFU)

Instructions for use are provided.

Conformity to European Regulations

This device is Category III PPE in accordance with regulation (EU) 2016/425 of the European Parliament and of the Council.

This device complies with the requirements of EN 166:2001 and associated standards and bear the following marks:

Marking	Description
1	Optical Class
3	Field of use: Liquid Product protects against liquid droplets or splashes
N	Resistance to fogging



Product Specification

Description	Pack Size	Box Size	Code
Reusable Innovia Face Visor	1x Visor	100x	CM-5001R

Material Breakdown

Product Component	Specification
Supportive Headrest	Adhesive Backed RH-32 Polyurethane Foam
Visor	250µm Polycarbonate, Clear
Headband	Latex Free Elastic
Individual Face Visor Packaging/ Bag	Grip Seal, Medium Square (12.75"x 12.75")
Outer Box/Carton Packaging	Plain Brown, Double Wall (610x457x457)



DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

DTR Medical Ltd, 17 Clarion Court, Enterprise Park, Swansea SA6 8RF, UK
Tel +44(0) 1792 79 79 10 Fax +44(0) 1792 79 79 55
marketing@dtmedical.co.uk | www.dtrmedical.com



NETWORK
MEDICAL PRODUCTS

AN
INNOVIA
MEDICAL
COMPANY

Coronet House, Kearsley Road, Ripon, North Yorkshire, HG4 2SG, UK
Tel: +44 (0) 1765 609555 Fax: +44 (0)1765 608476
info@networkmedical.co.uk | www.networkmedical.co.uk