



### Product Overview

The NETWORK ENT® Fine Ends range has been specifically designed to work alongside the NETWORK ENT® range of suction tubes to provide increased precision during delicate procedures. These bulk non-sterile products are designed for inclusion in custom procedure packs.

- Provides pin point suction for micro ENT procedures
- Manufactured with soft coloured-coded seal, for easy identification and superior fit
- Compatible with other makes of suction tube
- Supplied non-sterile, single use only
- Manufactured in the UK

### Fine End Size Options

Product	Size	Unit of Sale	Code
Fine End (Purple)	18G / 1.0 mm	Bag of 100	74-2020NS
Fine End (Orange)	20G / 0.8 mm	Bag of 100	74-2021NS
Fine End (Blue)	22G / 0.65 mm	Bag of 100	74-2022NS
Fine End (Green)	24G / 0.5 mm	Bag of 100	74-2023NS

### Material Specification

---

Product Component	Specification
Tube	Stainless Steel
Moulded handle	Polypropylene (PP)

---

### Intended Use

Network Aspirating and Suction tubes are designed particularly for ENT applications and when attached to a vacuum suction system, are used to remove excess fluids and debris from the surgical site. The Vari-Tube Fine ends are inserted into the tube ends of the Vari-Tube Suction handle to reduce the diameter of the tube for precision and pinpoint suctioning during delicate ontological procedures.

### Sterilisation

Products are supplied non-sterile, and are suitable for sterilisation by EO gas, gamma irradiation and e-beam. Devices are designated for Single Use Only.

### Instructions for Use

No instructions for use are provided.

### Conformity to the European Directives

The Network range of Aspiration Suction Tubes and Fine Ends are identified as surgically invasive devices for transient use and are therefore classified as a Class IIa devices according to Annex IX, Rule 6, of the Medical Devices Directive 93/42/EEC as amended by 2007/47/EC.