



PVA Eye Drains Specification Sheet PROCEDURE PACK (BULK) PRODUCT



Product Overview

The EYETEC® PVA Eye Drain range has been specifically designed for the quick and efficient management of fluids during ophthalmic procedures.

- Constructed from ultra-smooth, biocompatible PVA sponge
- Ultra-fast wicking
- Available sterile for direct inclusion in procedure packs
- Single use only, declared 5 year shelf life

PVA Eye Drain Options

Product	Pack Size	Unit of Sale	Product Code
Eye Drain 80cc Capacity	Foil of 1	Bag of 50 packs	2-430



PVA Eye Drains Specification Sheet PROCEDURE PACK (BULK) PRODUCT

Material Specification

Product Component	Specification
Eye Drain	100% Polyvinyl Alcohol Sponge (PVA), Poly bag
Pouch Packaging	Metalized foil pouch

Intended Use

The PVA Eye Drain range has been specifically designed for the quick and efficient management of fluids during ophthalmic procedures.

Intended Purpose

PVA eye drains:

A body and/or body orifice contact, single-use rectangular strip of medical grade, highly absorbent PVA which drains into an attached polythene bag which has a self-adhesive strip to permit it to be secured to the drape during ophthalmic surgical procedures. Intended to be used in a sterile condition by healthcare professionals on any patient for the management of fluids, and to staunch blood loss after invasive surgery or traumatic injury in the areas of ophthalmic surgery. They PVA portion of the drain has a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye), cornea and/or surrounding tissue. There are no known contraindications for these devices.



PVA Eye Drains Specification Sheet PROCEDURE PACK (BULK) PRODUCT

Sterilisation

Please note that PVA material is adversely affected by direct exposure to EO gas, therefore we supply these products inside a barrier foil that is sterilised by gamma irradiation. This allows the products to be included inside procedure packs that are sterilised by EO and maintains the optimum performance for the end-user.

Instructions for Use

The instructions for use are supplied in the form of a multi-lingual leaflet with instructions given. The international symbols used on the packaging are explained in each language. A leaflet is supplied with each product and is also available to view and download at <https://www.networkmedical.co.uk/ifu-product-group>

Conformity to the European Directives

EYETEC® PVA products are defined as invasive devices with respect to body orifices (Rule 5, Annex IX 93/42/EEC Medical Devices Directive) and are classified as Class I Sterile.